A.MENARINI
diagnostics

DETECTION DES AUTOANTICORPS

IVD

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REF 37783 Substrat Rein/estomac de souris 48 Tests

REF 37784 Substrat Rein/estomac/foie de souris 48 Tests

REF 37785 Lamelle Estomac/Rein de Souris 8 Puits

REF 37786 Lamelle Estomac/Rein/Foie de Souris 8 Puits

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REFERENCES • BIBΛIOΓPAΦIA • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Tan EM. Autoantibodies to nuclear antigens (ANA): Their immunobiology and medicine. *Adv Immunol* 33: 167-240, 1982.
2. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity of the skin. In "Concepts in Immuno-pathology", Vol 1, Cruse JM and Lewis RE Jr, Eds, Karger, Basel, 318-353, 1985.
3. Reimer G, Cornell RC and Tan EM. The biochemical nature of nuclear antigens reactive with antinuclear antibodies. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 519-531, 1987.
4. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Standardization of antinuclear antibody and other immunofluorescent tests used in immunopathologic studies of the skin. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 41-64, 1987.
5. Tan EM, Chan EKL, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): Diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. *Clin Immunol Immunopathol* 47: 121-141, 1988.
6. Manns M, Gerken G, Kyriatsoulis A and Meyer zum Büschenfelde KH. Significant autoimmune markers of autoimmune liver disorders: Current status. *J Clin Lab Anal* 1: 362-370, 1987.
7. Mackay IR. Autoimmunity and the liver. *Clin Aspects Immunity* 2: 8- 17, 1988.
8. McMillan SA, Alderdice JM, McKee CM et al. Diversity of autoantibodies in patients with anti-mitochondrial antibody and their diagnostic value. *J Clin Path* 4: 232-236, 1987.
9. Gershwin ME, Coppel RL and Mackay IR. Primary biliary cirrhosis and mitochondrial autoantigens - insights from molecular biology. *Hepatology* 8: 147-151, 1988.
10. Berg PA and Klein R. Mitochondrial antigens and autoantibodies from anti-M1 to anti-M9. *Klin Wochenschr* 64: 897-909, 1986.
11. Popper H and Paronetto F. Clinical, histologic and immunopathologic features of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 339-354, 1980.
12. Berg PA and Bacon H. Serology of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 355-373, 1980.
13. Anderson P, Small JV and Sobieszek A. Studies on the specificity of smooth muscle antibodies. *Clin Exp Immunol* 22: 22-29, 1975.
14. Kurki P, Miettinen A, Linder E, Pikkarainen P, Vuaristo M and Salaspuro MP. Different types of smooth muscle antibodies in chronic active hepatitis and primary biliary cirrhosis: Their diagnostic and prognostic significance. *Gut* 21: 878-884, 1980.
15. Fisher JB and Taylor KB. The significance of gastric antibodies. *Brit J Haematol* 20: 1-7, 1971.
16. Chisholm M. Immunology of gastritis. *Clin Gastroenterol* 5: 41 9-428, 1976.
17. Bigazzi PE, Burek CL and Rose NR. Antibodies to tissue-specific endocrine, gastrointestinal and neurological antigens. In "Manual of Clinical Laboratory Immunology". Rose NR, Friedman H and Fahey JL, Eds, American Society for Microbiology, Washington DC, 762-770, 1986.
18. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
19. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].
20. Nisengard RJ. Antinuclear antibodies: Significance of titers. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Bean S, Eds, John Wiley and Sons, New York, 2nd Ed, 387-398, 1979.

21. Meyer zum Büschenfelde KH, Manns M and Trautman F. Autoimmunity in chronic liver diseases - relationship to SLE? In "Recent Advances in Systemic Lupus Erythematosus". Lambert PH, Perrin L, and Izui S, Academic Press, New York, 259-269, 1984.
22. Walker JG, Doniach D, Roitt IM and Sherlock S. Serologic tests in the diagnosis of primary biliary cirrhosis. Lancet 1: 827, 1965.
23. Paronetto F and Popper H. Hetero-iso- and autoimmune phenomena in the liver. In "Textbook of Immunopathology", Miescher PA and Müller- Eberhard HJ, Eds, Grune and Stratton, New York, 2nd Ed, 789-817, 1976.
24. Leung PSC, Manns MP, Coppel RL, Gershwin ME. Detection of anti-mitochondrial antibodies in primary biliary cirrhosis and liver-kidney microsomal antibodies in autoimmune hepatitis. In "Manual of Clinical Laboratory Immunology", Rose NR, Hamilton RG, Detrick B, Eds, ASM Press, Washington DC, 6th Ed, 1023-1031, 2002.
25. Muratori P et al. Smooth muscle antibodies and type 1 autoimmune hepatitis. Autoimmunity. 35 (8): pp. 497–500. 2002.
26. Gatselis NK et al. Autoantibodies in HCV-treated patients. World J Gastroenterol. 11(4):482-487.2005.
27. Miller MH et al. Clinical comparison of cultured human epithelial cells and rat liver as substrates for the fluorescent antinuclear antibody test. J Rheumatol. 12 (2): 265-9. 1985.

Figure 1. Diagnostic. Significance of Antinuclear Antibodies

IF Staining Pattern	Nature of Antigen	Associated Disease
Homogeneous	Deoxyribonucleoprotein	SLE with renal involvement
Peripheral	DNA	SLE
Speckled	RNP	SLE or MCTD*
	Sm	SLE
	SS-A/SS-B	SLE or Sjögren's
	Sci-70	Scleroderma
Nucleolar	4S-6S RNA probably U3 RNA	Scleroderma
	Centromere/Kinetochores	inner and outer plates of kinetochores

*Mixed Connective Tissue Disease

Table 1: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on HEp-2 Cells

Clinical Condition	No. of Sera	% Positive
SLE	12	100
Subacute Cutaneous LE (SCLE)	7	86
Scleroderma	6	100
Rheumatoid Arthritis	10	50
Normal Controls	15	0

Table 2: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on Mouse Kidney Substrate

Clinical Condition	No. of Sera	% Positive
SLE	21	95
Scleroderma	17	82
Rheumatoid Arthritis	20	5
Normal Controls	96	0

Table 3: Incidence of Anti-Mitochondrial Antibodies (AMA) Detected by IFA on Mouse Kidney Substrate

Clinical Condition	% Positive
Primary Biliary Cirrhosis	100
Autoimmune Chronic Active Hepatitis	8
HBsAg and Chronic Active Hepatitis	0
Extrahepatic Jaundice and Other Liver Diseases	0
Systemic Lupus Erythematosus (SLE)	3
Rheumatoid Arthritis	0
Normal Controls	0

Adapted from Meyer zum Büschenfelde KH, et al.²¹; Walker JG, et al.²² and Paronetto F and Popper H²³.

Table 4: Incidence of Anti-Smooth Muscle Antibodies (ASMA) as Detected by IFA on Mouse Kidney Substrate

Clinical Condition	% Positive
Chronic Active Hepatitis (Type A)	50-87
Primary Biliary Cirrhosis	25
Acute Viral Hepatitis	87
Infectious Mononucleosis	87
Burkitt's Lymphoma	73
Nasopharyngeal Carcinoma	23
Hodgkin's Disease	23
Myeloproliferative Disorder	5
Warts	4
Normal Controls	3-18

Adapted from Anderson P, et al.¹³

Table 5. Incidence of Anti-Gastric Parietal Cell Antibodies (AGPA) as Detected by Indirect Immunofluorescence on Mouse Stomach Substrate

Clinical Condition	% Positive
Pernicious Anemia (PA)	85-95
Chronic Atrophic Gastritis without PA	30-60
Gastric Ulcer	25-30
Autoimmune Endocrinopathies	25-33
Sjögren's Syndrome	30
First Degree Relatives of PA Patients	30
Normal Controls	
< 20 years old	2
20-60 years old	6-8
> 60 years old	16

Table 6. Characteristics of LKM Antibody Positive Patients

	Condition	% Incidence
Sex (F/M)	4:1	
Age	Mean 34 years	
HBsAg	Negative	
Other Autoantibodies	Anti-nuclear	12%
	Anti-smooth muscle	18%
	Anti-mitochondrial	0%
	Anti-thyroglobulin	12%
	Anti-microsomal	35%
Immunoglobulin levels (gm/dl)	IgG	Normal
	IgM	Normal
	IgA	Low

Table 7. Findings in Positive Sera

	n	Negatives		Positives Titer	
		< 40	40-80	160-320	640-2560
ANA Positive Sera					
Menarini™	20	0	7	6	7
Other	20	0	7	8	5
AMA Positive Sera					
Menarini™	19	5	3	1	10
Other	19	5	3	4	7
ASMA Positive Sera					
Menarini™	19	7	8	4	0
Other	19	7	6	5	1
AGPA Positive Sera					
Menarini™	20	1	6	4	9
Other	20	2	6	7	5

Table 8. Findings in Normal Controls

	n	Negatives		Positives Titer	
		< 40	40-80	>160	
ANA Positive Sera					
Menarini™	38	33	5	0	
Other	38	36	2	0	
AMA Positive Sera					
Menarini™	38	38	0	0	
Other	38	38	0	0	
ASMA Positive Sera					
Menarini™	38	37	1	0	
Other	38	37	1	0	
AGPA Positive Sera					
Menarini™	38	37	1	0	
Other	38	37	1	0	

Table 9. Comparison of Kits Using Tissue Sections or HEp-2 Cells for the Detection of ANA

Clinical Condition	n	% Positive		
		Menarini™ Mouse Liver	Other Mouse Kidney	Other HEp-2
SLE	21	95	95	86
Scleroderma	17	82	88	77
RA	20	25	25	5
Normal Controls	96	0	0	0

A. Menarini Diagnostics S.r.l.
via Sette Santi 3
50131 Firenze
Italia

UK

UNITED KINGDOM

Distributed by

A. Menarini Diagnostics Ltd
405 Wharfedale Road
Winnersh - Wokingham
Berkshire RG41 5RA

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