# Test anticorps anti-ADN natif (ADNn)

Substrat *Crithidia luciliae*

IVD

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REF 37811 nDNA Kit 48 Tests

REF 37791 Lamelle Anti-nDNA *Crithidia luciliae* 6 Puits

## USAGE PRÉVU

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
## REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Tan EM. Autoantibodies to nuclear antigens (ANA): Their immunobiology and medicine. *Adv Immunol* 33: 167-240, 1982.
2. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity of the skin. In "Concepts in Immunopathology", Vol 1, Cruse JM and Lewis RE Jr, Eds, Karger, Basel, 318-353, 1985.
3. Reimer G, Cornell RC and Tan EM. The biochemical nature of nuclear antigens reactive with antinuclear antibodies. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 519-531, 1987.
4. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Standardization of antinuclear antibody and other immunofluorescent tests used in immunopathologic studies of the skin. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 41-64, 1987.
5. Tan EM, Chan EKL, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): Diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. *Clin Immunol Immunopathol* 47: 121-141, 1988.
6. Manns M, Gerken G, Kyriatsoulis A and Meyer zum Büschenfelde KH. Significant autoimmune markers of autoimmune liver disorders: Current status. *J Clin Lab Anal* 1: 362-370, 1987.
7. Mackay IR. Autoimmunity and the liver. *Clin Aspects Immunity* 2: 8- 17, 1988.
8. McMillan SA, Alderdice JM, McKee CM et al. Diversity of autoantibodies in patients with anti-mitochondrial antibody and their diagnostic value. *J Clin Path* 4: 232-236, 1987.
9. Gershwin ME, Coppel RL and Mackay IR. Primary biliary cirrhosis and mitochondrial autoantigens - insights from molecular biology. *Hepatology* 8: 147-151, 1988.
10. Berg PA and Klein R. Mitochondrial antigens and autoantibodies from anti-M1 to anti-M9. *Klin Wochenschr* 64: 897-909, 1986.
11. Popper H and Paronetto F. Clinical, histologic and immunopathologic features of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 339-354, 1980.
12. Berg PA and Bacon H. Serology of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 355-373, 1980.
13. Anderson P, Small JV and Sobieszek A. Studies on the specificity of smooth muscle antibodies. *Clin Exp Immunol* 22: 22-29, 1975.
14. Kurki P, Miettinen A, Linder E, Pikkarainen P, Vuaristo M and Salaspuro MP. Different types of smooth muscle antibodies in chronic active hepatitis and primary biliary cirrhosis: Their diagnostic and prognostic significance. *Gut* 21: 878-884, 1980.
15. Fisher JB and Taylor KB. The significance of gastric antibodies. *Brit J Haematol* 20: 1-7, 1971.
16. Chisholm M. Immunology of gastritis. *Clin Gastroenterol* 5: 419-428, 1976.
17. Bigazzi PE, Burek CL and Rose NR. Antibodies to tissue-specific endocrine, gastrointestinal and neurological antigens. In "Manual of Clinical Laboratory Immunology". Rose NR, Friedman H and Fahey JL, Eds, American Society for Microbiology, Washington DC, 762-770, 1986.
18. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
19. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].
20. Nisengard RJ. Antinuclear antibodies: Significance of titers. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Bean S, Eds, John Wiley and Sons, New York, 2nd Ed, 387-398, 1979.

21. Meyer zum Büschenfelde KH, Manns M and Trautman F. Autoimmunity in chronic liver diseases - relationship to SLE? In "Recent Advances in Systemic Lupus Erythematosus". Lambert PH, Perrin L, and Izui S, Academic Press, New York, 259-269, 1984.
22. Walker JG, Doniach D, Roitt IM and Sherlock S. Serologic tests in the diagnosis of primary biliary cirrhosis. Lancet 1: 827, 1965.
23. Paronetto F and Popper H. Hetero-iso- and autoimmune phenomena in the liver. In "Textbook of Immunopathology", Miescher PA and Müller-Eberhard HJ, Eds, Grune and Stratton, New York, 2nd Ed, 789-817, 1976.

**Figure 1. *C. luciliae***

**Table 1. Incidence of Antibodies to *n*DNA in Collagen-Vascular Disorders as detected by Indirect Immunofluorescence on *C. luciliae***

| Clinical Condition   | No. Tested | No. Positive | % Positive |
|----------------------|------------|--------------|------------|
| SLE                  | 28         | 19           | 68         |
| Scleroderma          | 23         | 0            | 0          |
| Rheumatoid Arthritis | 8          | 0            | 0          |
| Normal Controls      | 106        | 0            | 0          |

**Table 2. Comparison of Kits using *C. luciliae* Substrate for the Detection of Antibodies to *n*DNA**

| Clinical Condition   | n   | Menarini™ |            | Other    |            |
|----------------------|-----|-----------|------------|----------|------------|
|                      |     | Positive  | % Positive | Positive | % Positive |
| SLE                  | 28  | 19        | 68         | 13       | 46         |
| Scleroderma          | 23  | 0         | 0          | 0        | 0          |
| Rheumatoid Arthritis | 8   | 0         | 0          | 0        | 0          |
| Normal Controls      | 106 | 0         | 0          | 0        | 0          |

 **A. Menarini Diagnostics S.r.l.**  
via Sette Santi 3  
50131 Firenze  
Italia

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Benelux S.A./N.V.