Anti-human Tissue Transglutaminase (hu tTG) Screen IgG/IgA ELISA

IVD

PRODUCT INSERT

REF 38066 tTG Screen ELISA 96 Determinations

INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the detection and semiquantitation of anti-human Tissue Transglutaminase antibodies in human serum to aid in the diagnosis of celiac disease (CD) and dermatitis herpetiformis (DH).

SUMMARY AND EXPLANATION

Celiac Disease and Dermatitis Herpetiformis are autoimmune gastrointestinal disorders that may occur in genetically susceptible individuals triggered by the ingestion of gluten-containing grains such as wheat, barley and rye. CD is characterized by malabsorption resulting from inflammatory injury to the small intestinal mucosa and, when prolonged, can cause malnutrition. The classical symptoms of CD include diarrhea, weight loss and malnutrition. Only a small percentage of patients with CD, however, presents with classical symptoms. Studies have found the prevalence of CD to be highly variable from population to population¹. The true prevalence has been difficult to ascertain. If only the clinical criteria are used in determining prevalence, the incidence of CD is much lower as compared with incidence established by serological methods^{1,2}. Using serological methods the incidence of CD in the general population is approximately one in

200. Failure to diagnose CD early may predispose an individual to long-term complications such as splenic atrophy and intestinal lymphoma. A gluten-free diet (GFD) normalizes the mucosa and helps reduce the malignant potential^{3,4}. Histological examination of the small intestinal biopsy remains the gold standard for diagnosing CD, but it has limitations as well. For example, certain patients with latent or even active CD may have normal histopathology⁵. DH is differentiated from CD in that DH patients present with pruritic, itchy, vesiculo-bullous disorders of the skin in addition to underlying gluten enteropathy.

The revised European Society of Pediatric Gastroenterology and Nutrition (ESPGAN) criteria include only a single biopsy with clear-cut remission of clinical symptoms on GFD⁶. Positive serology at the time of diagnosis with disappearance on GFD contributes to the diagnosis. The various serological tests employed in the work-up of patients suspected to have CD include anti-gliadin antibody (AGA), anti-endomysial antibody (EMA), anti-reticulin antibody (ARA) and anti-tissue transglutaminase (tTG) antibody tests.

Since identification of tTG as the endomysial antigen, ELISA methods have been described for detecting antibodies in the sera of patients with CD. The advantage of the anti-tTG antibody assay is that it is automatable and less subjective than EMA. For this reason, many laboratories have opted to use the tTG antibody method as the screening assay. In various studies on the efficacy of the tTG antibody method for screening for CD, the specificity and sensitivity of this method has been found to range from 90 to 95%⁷. Human tTG has been described to improve the sensitivity of the tTG antibody assay for CD.

Most of the assays detect antibodies to tTG of the IgA isotype. This has some limitations as IgA deficient CD patients may yield false negative serology⁸. This may compromise the utility of the serum antibody methods in detecting all CD cases^{9,10}. IgA deficiency is a common immunodeficiency, found in one in 500-700¹¹⁻¹³

healthy blood donors. To prevent false negative results it is necessary to have a method that can detect antibodies to both immunoglobulin isotypes (IgA and IgG) against tTG.

PRINCIPLES OF PROCEDURES

The anti-hu tTG Total antibody test is performed as a solid phase immunoassay (ELISA). Microwells are coated with recombinant hu tTG antigen followed by blocking the unreacted sites to reduce nonspecific binding. Controls, calibrators and patient serum samples are incubated in the antigen coated wells, which allow anti-tTG antibodies present in the serum to bind. Unbound antibody and other serum proteins are removed by washing the microwells. Antibodies bound to the microwells are detected by adding enzyme labeled polyvalent anti-human immunoglobulin

conjugate to the wells. These enzyme conjugated antibodies bind specifically to the human immunoglobulin that are bound to tTG on the microwells. Unbound enzyme conjugate is removed by washing. Specific enzyme substrate (pNPP) is then added to the wells and the presence of antibodies to hu tTG is detected by a color change produced by the conversion of the pNPP substrate. The reaction is stopped and the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in ELISA units per milliliter (EU/ml).

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

Precautions

For *in vitro* diagnostic use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA licensed tests. However, human blood derivatives and patient specimens should always be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁴.

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources other than the same lot number from A. Menarini Diagnostics S.r.l.. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use beyond expiration date on the label.

Materials Provided

Menarini™ tTG Screen ELISA **REF** 38066

Kits contain sufficient reagents to perform 96 determinations each.

12 x 8	MICROPLATE hu-tTG	Microplate with individual breakaway microwells coated with hu tTG antigen.
1 x 1.5 ml	CALIBRATOR A hu-tTG Screen *	Ready to use Calibrator A (green cap). Human serum containing antibodies to tTG.
1 x 1.5 ml	CALIBRATOR B hu-tTG Screen *	Ready to use Calibrator B (violet cap). Human serum containing antibodies to tTG.
1 x 1.5 ml	CALIBRATOR C hu-tTG Screen *	Ready to use Calibrator C (blue cap). Human serum containing antibodies to tTG.
1 x 1.5 ml	CALIBRATOR D hu-tTG Screen *	Ready to use Calibrator D (yellow cap). Human serum containing antibodies to tTG.
1 x 1.5 ml	CONTROL + hu- tTG Screen *	Ready to use Positive Control (red cap). Contains human serum positive for tTG.
1 x 1.5 ml	CONTROL - *	Ready to use Negative Control (white cap). Contains human serum.

1 x 12 ml	IgA/IgG-CONJ	ALKPHOS *
1 x 60 ml	DIL *	
1 x 12 ml	SUBSTRATE *	
1 x 12 ml	STOP	
2 x	BUF	WASH

Ready to use **anti-human Alk. Phos. Conjugate**. Color coded pink.

Ready to use **Serum Diluent**. Color coded blue.

Ready to use **Enzyme Substrate**. Contains pNPP. **Protect from light**.

Ready to use **Stop Solution**.

Powder **Wash Buffer**. Reconstitute to one liter each.

* Contains < 0.1% NaN₃

Symbols used on labels:

LOT Lot number

REF Catalog number

Use by

Storage temperature

Read instructions for use

IVD In vitro diagnostic use

Manufacturer

Number of Tests

*Contains <0.1% NaN₃

Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettors capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Automatic microplate washer capable of dispensing 200 µl

Other available serological tests to aid in the Diagnosis of CD

Menarini™ Anti-Endomysial Antibody Test (Primate Distal Esophagus) (6x8 well slide) **REF** 38060

Menarini™ Anti-Endomysial Antibody Test (Primate smooth muscle) (8x6 well slide) **REF** 37790

Menarini™ Anti-gliadin Antibody (AGA) IgA ELISA **REF** 37789

Menarini™ Anti-gliadin Antibody (AGA) IgG ELISA **REF** 37809

Menarini™ Anti-gliadin Antibody (AGA) screen ELISA **REF** 37979

Menarini™ Anti-hu tTG Antibody IgA ELISA **REF** 37795

Menarini™ Anti-hu tTG Antibody IgG ELISA **REF** 37796

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freeze and thaw of samples.


PROCEDURE PROCEDURAL NOTES

- Before starting with the assay read carefully the product insert.
- Let serum specimens and test reagents equilibrate at room temperature for at least 30 minutes before starting with the test procedure. Return all unused specimens and reagents to refrigerator immediately after use.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Good washing technique is critical. If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.


Test Method

- Step 1** Let all reagents and specimens equilibrate at room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use only the Ready to Use Low Calibrator D (*vial with yellow cap*).
 or For a **semi-quantitative determination** use the Ready to Use Calibrators A through D as depicted in the sample layout below.

QUALITATIVE DETERMINATION

SEMI-QUANTITATIVE DETERMINATION

Specimen Layout

- Step 4** Prepare a **1:51** dilution of the patient samples by mixing **10 µl** of the patient sera with **0.5 ml** of Serum Diluent.
- Step 5** Pipette **100 µl** of Ready to Use Calibrators, Positive and Negative controls and diluted patient samples to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of the reagent blank should not be more than 0.3 when read against air.
- Step 6** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 7** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 8** Pipette **100 µl** of Conjugate into microwells.
- Step 9** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 10** Wash all microwells as in Step 7.
- Step 11** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 12** Incubate **microwells 30 minutes** (± 5 min) at room temperature.
- Step 13** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance within 1 hour of adding Stop Solution.
- Step 14** Read absorbance of each microwell at **405 nm** using a single or 405/ 630nm dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3 . Calibrator A should have an absorbance reading ≥ 1.0 , otherwise the test must be repeated. The negative control must be <20 EU/ml. If the test is run in duplicate, take the mean of the two readings to determine the concentration of anti-hu tTG. While performing qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial. We recommend borderline samples be tested with a fresh sample taken at a later date to ensure accuracy.

RESULTS

Calculations

Concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

Results obtained by this method should be reported as positive or negative.

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

2. SEMI-QUANTITATIVE DETERMINATION


Plot absorbance of Calibrator A through D against their respective concentration on linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of patient samples on the curve against the corresponding absorbance value.

Calibrator

The ready to use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such serum samples should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor.

Anti-hu tTG Total IgA/IgG Menarini™ Standard Curve

Interpretation

The following serves only as a guide in the interpretation of the laboratory results. These values were determined by testing 113 adult normal blood donors. The values depicted below are the mean of the normal subjects plus 3SD. Each laboratory must determine its own normal values.

anti-hu tTG value	Interpretation
<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

LIMITATIONS OF THE PROCEDURE

The Menarini™ anti-hu tTG assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. The method should be used for testing human serum samples only. Results obtained serve only as an aid in the diagnosis and should not be interpreted as diagnostic in themselves.

EXPECTED VALUES

The expected values in a normal population are negative (<20 EU/ml). In patients with CD the incidence of IgA anti-tTG antibodies has been described in 82-100% of patients with active CD with a specificity of 76-99%¹⁵. Combined detection of IgG and IgA anti-tTG antibodies as described by Seissler et al provided 95% sensitivity with 99% specificity¹⁶. The incidence and the levels of anti-tTG antibodies is dependent upon the diet status. The levels of these antibodies decrease and eventually will become negative in patients with CD who are on a gluten-free diet. Similarly, the levels of these antibodies will increase and may become positive when patients with CD who were on a gluten-free diet ingest a gluten-containing diet¹⁷⁻²¹.

PERFORMANCE CHARACTERISTICS

The utility of the Menarini™ anti-hu tTG Total Antibody ELISA was determined by comparing the results with:

- Clinical findings.
- Comparison to endomysial antibody (EMA)

Normal Range: The normal range was established by testing 65 serum samples from apparently healthy donors obtained from the Red Cross. The mean plus three standard deviations of the mean of this normal population was used to determine the cut-off between normal and borderline positive individuals.

Comparative Specificity and Sensitivity:

A. Clinical Correlation: A total of 200 serum specimens obtained on patients with CD with and without IgA deficiency; various disease controls including sera from patients with clinical diagnosis of the blistering diseases pemphigus, pemphigoid, and another dermatological disorder, psoriasis; and normal subjects were tested for total anti-hu tTG antibodies. The incidence of total anti-tTG antibodies in CD and other disease controls is as below. Four normals who were found positive were weakly positive.

Diagnosis	No Tested	No. Positive
Celiac Disease	94	75*
IgA deficient		
Celiac	15	15
Non-Celiac	13	0
Dermatitis herpetiformis	22	14
Disease Controls	56	2
Normal	65	4

*The 19 Celiac Samples that were negative for total tTG, were tested on both tTG IgG ELISA and tTG IgA ELISA. In all cases the 19 samples were negative on the two different ELISA tests. Most likely, these patients were on a gluten free diet.

B. Total anti-hu tTG ELISA vs. -endomysial antibody (EMA): A total of 200 samples were tested on the Menarini™ total anti-hu tTG antibody kit and EMA test kit. Two of the four samples positive on Menarini™ total tTG and negative on the EMA were sera from patients with DH and were weak positives. The other two false positives were within the Pemphigus/Pemphigoid disease population set. Three of the five samples negative on Menarini™ total tTG kit and positive on EMA had titers of 10 or less. The other samples in the CD subset had titers of 40 and 20 with 14.6EU/ml and 19.8EU/ml results respectively.

		Menarini™ EMA		
		Positive	Negative	Total
Menarini™ anti-hu total tTG	Positive	102	4	106
	Negative	5	89	94
	Total	107	93	200

Relative Agreement: 96%

Relative Sensitivity: 95%

Relative Specificity: 96%

C. Cross Reactivity: A total of 56 controls (patients with connective tissue disorders such as psoriasis, autoimmune bullous disorders such as pemphigus and pemphigoid) were tested. Two were found weakly positive.

Precision:

Based on 10 replicates, the intra-assay and inter-assay Coefficient of Variation (CV) of the total anti-hu tTG ELISA test were calculated.

	Inter-assay CV	intra-assay CV
High	5.5%	3.8%
Medium	9.7%	5.8%
Low	8.5%	6.4%
Cut-off (20 EU/ml)	11%	9%


Linearity:

To determine acceptable linearity, plates were assayed with calibrators of known values. The r-squared values of the standard curve were determined. An r-squared value greater than 0.95 is deemed to be acceptable. The average r-squared value for this assay 0.9975. No value was lower than 0.9929.

Recovery:

Samples with known IgG anti-hu tTG antibody concentrations were mixed with appropriate dilutions of another positive sample with known amounts of anti-tTG antibody. Anti-hu tTG antibody levels of the mixed samples were determined and from the values obtained the percent recovery calculated. The results are as follows:

	Anti-hu tTG Abs. conc. added (EU/ml)	Anti-hu tTG Abs. conc. obtained (EU/ml)	% Recovery
Sample 1	54.1	53.7	99.4
Sample 2	73.8	86.5	117.2
Sample 3	89.3	106.1	118.7

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
RESUMO DO PROCEDIMENTO Menarini™

REFERENCES • BIBΛIOΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Guandalini S, Gupta P. Celiac disease- A diagnostic challenge. *Clin Appl Immun Rev* 2002;2:293-305.
2. Fasano A, et al. Prevalence of celiac disease in at risk and not-at-risk groups in the United States 2003; 163:286-292.
3. Wolber R, et al. Lymphocytic gastritis in patients with celiac sprue or spruelike intestinal disease. *Gastroenterology* 1990;98:310-315.
4. Catassi, et al. Risk of non-Hodgkin lymphoma in celiac disease. *JAMA* 2002;287:1413-1419.
5. Arranz E, Ferguson A. Intestinal antibody pattern of celiac disease: Occurrence in patients with normal jejunal biopsy histology. *Gastroenterology* 1993;104:1263-72.
6. Hill ID, et al. Celiac disease: Working group report of the first world congress of pediatric gastroenterology, hepatology, and nutrition. *J. Pediatr Gastroenterol Nutr* 2002; 35:S78-S88
7. Vitoria JC, et al. Antibodies to gliadin, endomysium, and tissue transglutaminase for the diagnosis of celiac disease. *J Pediatr Gastroenterol Nutr* 1999;29:571-4.
8. Rittmeyer C, Rhoads JM. IgA deficiency causes false-negative endomysial antibody results in celiac disease. *J Pediatr Gastroenterol Nutr*. 1996 23:504-6.
9. Kumar V, Jarzabek-Chorzelska M, Sulej J, Rajadhyaksha M, Jablonska S. Tissue transglutaminase and endomysial antibodies-diagnostic markers of gluten-sensitive enteropathy in dermatitis herpetiformis. *Clin Immunol*. 2001;98:378-82.
10. Lagerqvist, C. Ivarsson, A., Juto, P. Persson, L. A. Hernell, O. Screening for adult coeliac disease - which serological marker(s) to use? *J Intern Med* 250:241-248, 2001.
11. Lilic D, Sewell WA. IgA deficiency: what we should-or should not-be doing. *J Clin Pathol* 2001 54:337-8.
12. Seidman EG, Hollander GA. Autoimmunity with immunodeficiency: a logical paradox. *J Pediatr Gastroenterol Nutr*. 1999 ;28:377-9.
13. Sleasman JW. The association between immunodeficiency and the development of autoimmune disease. *Adv Dent Res* 1996;10:57-61.
14. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).
15. Korponay-Szabo IR, Dahlbom I, Laurila K, et al. Elevation of IgG antibodies against tissue transglutaminase as a diagnostic tool for coeliac disease in selective IgA deficiency. *Gut*. 2003;52:1567-1571.
16. Leon F, R-Pena C, Sanches L et al. Anti-transglutaminase IgA ELISA: Clinical potential and drawbacks in celiac disease diagnosis. *Scand J Gastroenterol* 2001;36:849 853.
17. Seissler J, Boms S, Wohlrab U et al. Antibodies to human recombinant tissue transglutaminase measured by radioligand assay: evidence for high diagnostic sensitivity for celiac disease. *Horm Metab Res* 1999;31:375-379.
18. Kumar V et al. Celiac disease and immunoglobulin a deficiency: how effective are the serological methods of diagnosis? *Clin Diagn Lab Immunol*. 2002; 9:1295-300.
19. Cataldo F, Lio D, Marino V, Picarelli A, Ventura A, Corazza GR. IgG antiendomysium and IgG antitissue transglutaminase (anti-tTG) antibodies in coeliac patients with selective IgA deficiency. Working Groups on Celiac Disease of SIGEP and Club del Tenue. *Gut*. 2000 47:366-9.
20. Sulkanen S, Haltunen T, Laurila K, et al. Tissue Transglutaminase autoantibody enzyme-linked immunosorbent assay in detecting celiac disease. *Gastroenterology*; 1998, 115:1322-1328.
21. Dietrich W, Laag E, Schopper H, et al. Autoantibodies to tissue transglutaminase as predictors of celiac disease; *Gastroenterology* 1998, 115:1317-1321.

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Benelux S.A./N.V.