# Oxidized LDL (oxLDL) Antibody ELISA

IVD

## PRODUCT INSERT

REF 38085 Oxidized LDL Antibody (oxLDL) ELISA 96 Determinations

### INTENDED USE

Enzyme linked immunosorbent assay (ELISA) for the qualitative and semi-quantitative detection of antibodies to oxidized low-density lipoprotein (oxLDL) in human serum.

### SUMMARY AND EXPLANATION

Anti-oxidized LDL antibodies are associated with certain clinical conditions such as atherosclerosis, carotid stenosis and other coronary artery dysfunction. Anti-oxidized LDL antibodies may also be significant in diseases such as diabetes hypertension and endometriosis.<sup>1-5</sup>

Hypercholesterolemia is closely associated with increased risk of atherosclerosis.<sup>4,6</sup> Cholesterol accumulated in the atherosclerotic plaque is derived primarily from low-density lipoprotein (LDL). Oxidation of LDL is accepted as a critical event in the atherogenic process.<sup>7</sup> Free radicals like superoxide and nitric oxide ( $-O_2$ , NO) generated in biological reactions in the body contribute to the oxidation of LDL.<sup>8,9</sup> The NO radical has been shown to oxidize apolipoprotein B, a constituent of LDL. Similarly, lipo-oxygenases and oxidants like per-oxynitrite can oxidize the lipid moieties in LDL. It has been demonstrated that oxidized LDL and not native LDL is taken up by scavenger receptors on monocytes, smooth muscle cells and macrophages in the blood vessels. The oxidation of LDL increases the affinity of oxidized LDL to acetyl receptors of macrophages.<sup>10</sup> As this pathway of oxidized LDL uptake is unregulated, it results in the formation of lipid-laden macrophages (foam cells). The presence of foam cells is the earliest step in the formation of atherosclerotic plaque in blood vessels. Later, other inflammatory molecules released by leukocytes promote the progression of atherosclerotic plaque.<sup>11,12</sup>

Antibodies against oxLDL and its complexes have been demonstrated in systemic lupus erythematosus, diabetes, arterial thrombosis, hypertension, etc.

### PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with oxLDL antigen followed by a blocking step to reduce non-specific protein binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies in samples to bind to the antigen. Unbound antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgG conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in ELISA units per milliliter (EU/ml) and reported as positive or negative.

### REAGENTS

#### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.**

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Coated microwell strips are for one time use only.

## Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.<sup>13</sup>

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

## Materials provided

Zenit oxLDL ELISA [REF] 38085

Kits contain sufficient reagents to perform 96 determinations.

12 x 8	MICROPLATE oxLDL	Microplate with individual breakaway microwells. Coated with oxidized LDL antigen. Ready for use.
1 x 1.75 ml	CONTROL + oxLDL	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for oxLDL antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL -	Ready to use <b>Negative Control</b> ( <i>white cap</i> ). Contains human serum.
5 x 1.75 ml	CALIBRATOR A oxLDL CALIBRATOR B oxLDL CALIBRATOR C oxLDL CALIBRATOR D oxLDL CALIBRATOR E oxLDL	Ready to use <b>set of 5 Calibrators</b> . Derived from human serum containing oxLDL antibodies. Concentrations in EU/ml are printed on the labels. Calibrator A ( <i>green cap</i> ) 320 EU/ml; Calibrator B ( <i>violet cap</i> ) 80 EU/ml; Calibrator C ( <i>blue cap</i> ) 40 EU/ml; Calibrator D ( <i>yellow cap</i> ) 20 EU/ml, and Calibrator E ( <i>orange cap</i> ) 1 EU/ml.
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human IgG. Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded green.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. <b>Protect from light.</b>
1 x 15 ml	STOP	Stop Solution. Ready for use.
2 x	BUF WASH	Powder Wash Buffer. <b>Reconstitute to one liter each.</b>

## Symbols used on labels

LOT	Lot number
REF	Catalog number
IVD	In vitro diagnostic use
	Use by
	Storage temperature
	Read instructions before use
	Number of tests
	Manufacturer

### Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

### PROCEDURE


#### Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.
- **Ensure microplate incubations occur in the dark. This includes specimen, conjugate and substrate incubation steps. Exposure to light through incubation may impact results.**


#### Test Method

- Step 1** Let all reagents and specimens equilibrate to room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use Calibrator D (*yellow cap*) only.  
or  
For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

## Qualitative

## Semi-quantitative

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500 µl** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.  
**Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** ( $\pm$  5 min) at room temperature in the dark.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** ( $\pm$  5 min) at room temperature in the dark.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** ( $\pm$  5 min) at room temperature in the dark in the dark.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within 1 hour of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be  $<0.3$ . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be  $<20$  EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

## RESULTS

### Calculations

The concentrations of the samples can be determined by either of two methods:

#### 1. QUALITATIVE DETERMINATION

Abs. of Test Sample

$$\text{-----} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

Abs. of Calibrator D

#### 2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

### Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

### Interpretation

Interpretation values were determined by testing 64 adult normal blood donors. The mean of the normal subjects plus 2 SD was established as the assay cutoff and assigned an arbitrary value of 20 EU/ml. A. Menarini Diagnostics suggests use of the reference range below. Each laboratory should validate assay values for their own conditions.

<b>oxLDL Ab value</b>	<b>Interpretation</b>
<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

### LIMITATIONS OF PROCEDURE

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

Results obtained by this assay alone are not diagnostic and should be considered in conjunction with clinical presentation of the patient. Any borderline patient result should be retested to confirm the result. It is also recommended that patients with borderline results be retested at appropriate intervals.

### EXPECTED VALUES

Test results in a normal population are expected to be negative. Oxidized LDL antibodies are present in patients with atherosclerosis, diabetes, SLE and stress inducible disorders.

### Performance Characteristics

A. Disease association: The Zenit oxLDL Antibody ELISA for the detection of oxidized LDL antibodies was evaluated by testing sera of patients with associated conditions as reported in the literature, disease control sera

from patients with unassociated autoimmune diseases and control sera from healthy normal subjects. Results of these tests and the number of sera taken from each population are summarized below.

Population	n	oxLDL Ab Pos n (%)
Coronary Artery Disease	9	7 (78%)
Phospholipid Antibody Positive	57	33 (58%)
Disease Controls	50	1 (2%)
Normals	100	4 (4%)

B. Cross Reactivity: A total of 50 potentially cross-reactive specimens from individuals with other autoimmune disorders or positive for other autoantibodies were tested for oxLDL antibodies using the Zenit oxLDL system.

Condition	n	oxLDL Ab positive n (%)
Antinuclear Ab positive	10	0 (0%)
Cyclic Citrullinated Peptide Ab positive	10	0 (0%)
Hashimoto's disease	10	0 (0%)
Hsp-70 Ab positive	10	0 (0%)
Rheumatoid Factor positive	10	1 (10%)
Total	50	1 (2%)

### Precision

Inter-assay precision was tested with 3 positive specimens selected throughout the range of the assay. Multiple assay runs were performed at separate times. Intra-assay precision was calculated based on 10 replicates of 3 specimens. Based on these replicates the Coefficient of Variation (CV) was calculated. Inter-assay and intra-assay CV ranged from 5-13%.

### Linearity

To determine acceptable linearity, plates were assayed with calibrators of known values. The r-squared values of the resulting standard curves were determined for multiple assays. An r-squared value greater than 0.95 is deemed to be acceptable. The average r-squared value for these assays in precision testing was 0.9977 with no value lower than 0.9972.

### Recovery

Samples with known oxLDL antibody levels were mixed with appropriate dilutions of other positive samples with known amounts of oxLDL antibodies. The antibody levels of the mixed samples were determined and used to calculate percent recovery. The results follow:

Samples	S1 (EU/ml)	S2 (EU/ml)	EU/ml Expected	EU/ml Obtained	% Recovery (Obtained/Expected)
Mixture 1	43.6	42.2	42.9	49.4	115.2%
Mixture 2	37.1	60.2	48.6	50.4	103.6%
Mixture 3	20.1	52.9	36.5	38.0	104.0%
Mixture 4	78.9	65.3	72.1	84.3	116.9%
Mixture 5	178.5	152.8	165.7	171.3	103.4%

### Interference

Interference was studied by mixing sera with known oxLDL antibody levels with potentially interfering serum samples and studying deviation from expected results. Results appear below.

Mixture*	S1 (EU/ml)	S2 (EU/ml)	EU/ml Expected	EU/ml Obtained	%Interference
H1	24.0	19.1	21.6	24.1	111.8%
L1	24.0	14.0	19.0	19.9	104.9%
H2	42.3	19.1	30.7	26.8	87.2%
L2	42.3	14.0	28.1	24.0	85.4%
H3	79.3	21.5	50.4	48.3	95.9%
L3	79.3	18.6	48.9	48.9	99.9%
H4	189.8	12.2	101.0	131.8	130.5%
L4	189.8	2.1	95.9	117.1	122.1%
H5	167.9	24.6	96.3	101.3	105.3%
L5	167.9	11.6	89.8	104.1	115.9%

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