Dr. Fenning
BioMed GmbH

fr-elisa

Autoimmune-Diagnostic

β 2-Glycoprotein 1-IgA

11024

Content:

English	2
Italian	4
German	6
Spanish	8
Portuguese	10
French	12
Greek	14
Technical Data	16
Explanation of symbols	17

Enzyme Immunoassay for the Detection of β₂-Glycoprotein 1-IgA

Product-No. 11024

Important Comments:

- Before use, please carefully read the instructions for the fr-elisa.
- The kit should be stored at 2-8°C. Do not freeze.
- Please notice the expiry date of all individual reagents and do not use them if overdue.
- Warm up the contents of the testkit to RT!
- Always open the foil pouch of the microtiterplate after reaching RT.
- Washing buffer and sample buffer are presented in a higher concentrated form compared to the working dilutions required in the test. After diluting the concentrated solutions, please regard the following use-by dates. Sample buffer as well as washing buffer diluted (stored at 2° - 8°C): 2 weeks.
- Do not mix up reagents of different fr-elisa Test Kit lots.

Clinical Significance:

New evidence is emerging that the presence of anti-phospholipid antibodies may not be the only cause involved in autoimmune thrombotic disorders. The binding of anti-cardiolipin antibodies was recently demonstrated to depend on the presence of cofactor β₂-glycoprotein 1 (β₂-GP 1), also known as apolipoprotein H. β₂-glycoprotein 1 is a 326-amino acid protein with an apparent molecular weight of 50 kD. Some studies suggest that β₂-glycoprotein 1 is the true antigen target for anti-phospholipid antibodies, whereas others postulate the formation of a cardiolipin / β₂-glycoprotein 1 complex as antigen target. Nevertheless, the binding of β₂-glycoprotein 1 to anionic lipids or surfaces is generally considered essential to differentiate between anti-cardiolipin antibodies associated with autoimmune diseases and anti-cardiolipin antibodies associated with infectious diseases like tuberculosis, leprosy, Lyme disease, AIDS, mononucleosis, hepatitis, etc. β₂-glycoprotein 1 without phospholipids is more specific to detect potential venous or arterial thrombosis. This specificity has been demonstrated by a close correlation between lupus anticoagulant and anti-β₂-glycoprotein 1 antibodies. In this study 88% of sera positive for anti-β₂-glycoprotein 1 antibodies also showed lupus anticoagulant activity, whereas only 18% of sera positive for lupus anticoagulant activity were negative for anti-β₂-glycoprotein 1 antibodies. The β₂-glycoprotein 1 autoantibody test is therefore more specific and it is useful to use it in addition to the traditional anti-cardiolipin and lupus anticoagulant tests for assessing the risk of thrombosis in a risk population.

- MC NEIL HP, SIMPSON RJ, CHESTERMAN CN and KRILIS SA: 1990 *Proc Natl Acad Sci, USA* 87: 4120-4124
- GALLI M, COMFUTIUS P, MAASSEN C, HEMKER HC, DE BAETS MH, VAN BREDA-VRIESMAN PJC, BARBUI T, ZWAAL RFA and BEVERS EM: 1990 *Lancet* 335: 1544-1547
- MATSUURA E, IGARASHI Y, FUJIMOTO M, ICHIKAWA K and KOIKE K: 1990 *Lancet* 336: 177-178
- ERICSON E, NAJMEY S, KEIL L, EL-KADI H and DE BARI: 1996 *Clinical Chemistry*, 42: 1-2
- VIARD JP, AMOURA Z and BACH JF: 1992 *Am J Med.* 93: 181-186

Principle of the Test:

The kits are solid-phase enzyme immunoassays. Microtiter wells are coated with affinity-purified human β₂-glycoprotein 1, cardiolipin cofactor. Antigen-precoated microplate wells are incubated with calibrators, controls and serum specimens. During the incubation, antibodies present in the test sample bind to the coated wells. Horseradish peroxidase-conjugated anti-human IgA is incubated in the wells to recognize the autoantibodies bound to the coated wells. At the end of each incubation, the unbound material is removed by aspirating and washing. Chromogen is added and autoantibodies are measured using a spectrophotometric plate reader.

Contents of the Test Kit:

- 12 x 8 coated microtiter test strips, prefixed in the frame
- 4 x calibrators, ready to use 3.125 SA U/ml; 12.5 SA U/ml; 50 SA U/ml; 100 SA U/ml; 1000µl each. **Note: The calibration curve has to be prepared for 4 + 1 calibrators; the (prediluted) sample dilution buffer is used as "zero calibrator".**
- 1 x negative control serum; ready to use; 1000µl
- 1 x positive control serum; ready to use; 1000µl
- 1 x sample dilution buffer (**ochre** coloured); 5-fold concentration; 20ml
- 1 x enzyme conjugate solution (**green** coloured); 15ml
- 1 x washing buffer; 20-fold concentration; 50ml
- 1 x substrate solution (TMB); 15ml
- 1 x stop solution, 10ml
- data sheet: measuring range of control sera of the test
- product information
- pipetting scheme

Test Performance:

We recommend the usage of multi-channel pipette and an automatic washer to achieve high synchronous incubation times and high performance for all calibrators and samples used in the assay.

Preliminary steps:

- operate the test kit at room temperature
- always open the foil pouch after reaching RT
- dilute the sample dilution buffer **ochre** (1 part) with distilled water (4 parts). Total volume = 100ml
- dilute the washing buffer (1part) with distilled water (19 parts). Total volume = 1000 ml
- the patient sera must be diluted 1:101 prior to use.
Dispense 10µl specimen into 1 ml diluted sample buffer **ochre**.
- the calibrators are ready to use
- the positive and negative control sera are ready to use
- the enzyme conjugate solution **green** is ready to use
- the chromogen (TMB) has to be used without any influence of light. After usage it has to be stored protected from light
- Unused test strips should be promptly resealed in the foil pouch with desiccant and stored at 2-8°C.
- **Handle the stop solution carefully! Sulphuric Acid**

Test Procedure:

- **100µl of the serum dilution**, the undiluted **calibrators** and the undiluted **control sera** are pipetted into the wells. It is recommended to perform repeat determinations and "blanks" (sample dilution buffer **ochre**, pre-diluted)
- Incubate **30 minutes** at room temperature
- rinse the microtiter plate three times using at least 300µl washing buffer for each well and each washing step. Remove traces of remaining buffer out of the wells by blotting thoroughly on absorbent paper after the last wash. Attention: incomplete washing and aspiration of wells may lead to decreased precision.
- add **100µl of enzyme conjugate** solution **green** to each well
- incubate **30 minutes** at room temperature
- discard enzyme tracer from wells and wash plate as described above.
- add **100µl of chromogen** (TMB) to each well
- incubate **30 minutes** completely protected from light at room temperature
- add **50µl of stop solution** in the same order as chromogen
- determine the optical extinction at **450nm** within 30 min after completing the assay.

Evaluation of the Test:

The extinctions determined for the calibration points are marked versus the concentrations of the standard solutions semilogarithmically (3.125 SA U/ml; 12.5 SA U/ml; 50 SA U/ml; 100 SA U/ml). A calibration curve is thus obtained by drawing a best-fit line (four parameters) through the points. Directly from the calibration curve, read the concentration of each analyte.

β2-Glycoprotein 1-IgA values lower than 4 SA U/ml are judged as negative.

Cut-off = 4 SA U/ml

Interpretation of results depends on the specific clinical application of the test: each laboratory should establish its own clinically relevant ranges for the population taken into consideration. A specimen with equivocal antibody levels should be tested again; if it remains equivocal, the result should be reported as equivocal and/or an additional sample should be taken for testing according to physician judgment.

The test can be evaluated if the positive control serum is determined in the range given by the data sheet and if in the same time the negative control is found below the "Cut-Off" value.

Precautions:

For in-vitro diagnostics only! Standard and control sera are of human origin. The sera were tested and found to be negative for HBsAg, Hepatitis C, and HIV. Nevertheless all reagents should be regarded as potential infectious and thus handled with the care required. Rules for handling human sera have to be obeyed.

Warning: Some reagents contain Sodium Azide. Sodium Azide can form explosive Metal Azides with plumb and copper. Residuals of reagents should be removed with water carefully. Some reagents contain small amounts of Thimerosal (<0.1% w/v). Substrate contains 3,3', 5,5' Tetramethylbenzidine. The stop solution contains 2,6% Sulphuric Acid. Therefore handle all components as if potentially hazardous.

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CE Dr. Fenning BioMed GmbH