Dr. Fenning
BioMed GmbH

fr-elisa

Autoimmune-Diagnostic

Cardiolipin-IgA

11027

Content:

English	2
Italian	4
German	6
Spanish	8
Portuguese	10
French	12
Greek	14
Technical Data	16
Explanation of symbols	17

Enzyme Immunoassay for the Detection of Cardiolipin-IgA

Product-No. 11027

Important Comments:

- Before use, please carefully read the instructions for the fr-elisa.
- The kit should be stored at 2-8°C. Do not freeze.
- Please notice the expiry date of all individual reagents and do not use them if overdue.
- Warm up the contents of the testkit to RT!
- Always open the foil pouch of the microtiterplate after reaching RT.
- Washing buffer and sample buffer are presented in a higher concentrated form compared to the working dilutions required in the test. After diluting the concentrated solutions, please regard the following use-by dates. Sample buffer as well as washing buffer diluted (stored at 2° - 8°C): 2 weeks.
- Do not mix up reagents of different fr-elisa Test Kit lots.

Clinical Significance:

Autoantibodies against phospholipids (APAs) have been associated with the occurrence of thrombosis, thrombocytopenia, fetal abort, cerebro-vascular insufficiency and myocardial infarction (anti-phospholipid syndrome). Anti-cardiolipin antibodies may be IgG, IgM or IgA: the clinical significance of the three main immunoglobulin subclasses is a subject of controversial discussion. The presence of IgG anti-cardiolipin antibodies appears to be the most predictive and specific test for thrombosis, thrombocytopenia, and recurrent fetal loss. In addition anti-cardiolipin IgG antibody titers are closely associated with neurologic complications in systemic lupus erythematosus patients and in patients with other autoimmune diseases. However patients with anticardiolipin-associated clinical features often have elevated IgM and IgA anti-cardiolipin antibodies. Therefore it has become clinically relevant to determine the presence of each of these classes of autoantibodies in patient serum.

- CRONIN ME, BISWAS RM, VAN DER STRAETON C, FLEISCHER TA et al.: IgG and IgM anticardiolipin antibodies in lupus anticardiolipin antibody associated clinical syndromes. J. Rheumatol., 15: 795-798 (1988).
- COWCHOCK S, FORT J, MUNOZ S, NORBERG R et al. False positive ELISA test for anticardiolipin antibodies in sera from patients with repeated abortions, rheumatologic disorders, and primary biliary cirrhosis: Correlation with elevated polyclonal IgM and implications for patients with repeated abortions. Clin. Exp. Immunol., 73: 289-294 (1988).
- HARRIS EN. Solid-phase anticardiolipin test revisited [editorial], Am. J. Med., 85: 599-601 (1988).
- WILSON WA; GHARAVI AE, KOIKE T, LOCKSHIN MD et al. International consensus statement on preliminary classification criteria for definite antiphospholipid syndrome. Arthritis & Rheumatism Vol.42 No.7: 1309-1311 (1999)
- LOCKSHIN MD, SAMMARITANO LR, SCHWARTZMAN S. Validation of the Sapporo criteria for antiphospholipid syndrome. Arthritis & Rheumatism Vol.43 No.2: 440-443 (2000)
- MIYAKIS S, LOCKSHIN MD, ATSUMI T, BRANCH DW, BREY RL, CERVERA R et al. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). J. Thromb and Haemost 4: 295-306 (2006)

Principle of the Test:

The kits are solid-phase enzyme immunoassays. Microtitre wells are coated with cardiolipin and β 2-glycoprotein I (β 2-GPI, apolipoprotein H), the latter being cardiolipin co-factor, which forms complexes with cardiolipin. The presence of cardiolipin and β 2-glycoprotein 1 complexes allows specific cardiolipin antibodies to bind to the solid phase. Antigen-precoated microplate wells are incubated with calibrators, controls and serum specimens. During the incubation, antibodies present in the test sample bind to the coated wells. Horseradish peroxidase-conjugated anti-human IgA is incubated in the wells to recognize the autoantibodies bound to the coated wells. At the end of each incubation the unbound material is removed by aspirating and washing. Chromogen is added and autoantibodies are measured using a spectrophotometric plate reader.

The calibrator concentrations are set up to: 3.125, 12.5, 50, 100 APL, referenced to Prof. N. Harris Reference Standard (Louisville University, KY, USA). 1 APL corresponds to 1 μ g/ml of affinity-purified IgA.

Contents of the Test Kit:

- 12 x 8 coated microtiter test strips, prefixed in the frame
- 4 x calibrators, ready to use 3.125APL; 12.5 APL; 50 APL; 100 APL; 1000 μ l each. **Note: The calibration curve has to be prepared for 4 + 1 calibrators; the (prediluted) sample dilution buffer is used as "zero calibrator".**
- 1 x negative control serum; ready to use; 1000 μ l
- 1 x positive control serum; ready to use; 1000 μ l
- 1 x sample dilution buffer (**ochre** coloured); 5-fold concentration; 20ml
- 1 x enzyme conjugate solution (**green** coloured); 15ml
- 1 x washing buffer; 20-fold concentration; 50ml
- 1 x substrate solution (TMB); 15ml
- 1 x stop solution, 10ml

- data sheet: measuring range of control sera of the test
- product information
- pipetting scheme

Test Performance:

We recommend the usage of multi-channel pipette and an automatic washer to achieve high synchronous incubation times and high performance for all calibrators and samples used in the assay.

Preliminary steps:

- operate the test kit at room temperature
- always open the foil pouch after reaching RT
- dilute the sample dilution buffer **ochre** (1 part) with distilled water (4 parts). Total volume = 100ml
- dilute the washing buffer (1part) with distilled water (19 parts). Total volume = 1000 ml
- the patient sera must be diluted 1:101 prior to use.
Dispense 10µl specimen into 1 ml diluted sample buffer **ochre**.
- the calibrators are ready to use
- the positive and negative control sera are ready to use
- the enzyme conjugate solution **green** is ready to use
- the chromogen (TMB) has to be used without any influence of light. After usage it has to be stored protected from light
- Unused test strips should be promptly resealed in the foil pouch with desiccant and stored at 2-8°C.
- **Handle the stop solution carefully! Sulphuric Acid**

Test Procedure:

- **100µl of the serum dilution**, the undiluted **calibrators** and the undiluted **control sera** are pipetted into the wells. It is recommended to perform repeat determinations and "blanks" (sample dilution buffer **ochre**, pre-diluted)
- **Incubate 30 minutes** at room temperature
- rinse the microtiter plate three times using at least 300µl washing buffer for each well and each washing step. Remove traces of remaining buffer out of the wells by blotting thoroughly on absorbent paper after the last wash. Attention: incomplete washing and aspiration of wells may lead to decreased precision.
- add **100µl of enzyme conjugate** solution **green** to each well
- **incubate 30 minutes** at room temperature
- discard enzyme tracer from wells and wash plate as described above.
- add **100µl of chromogen** (TMB) to each well
- **incubate 30 minutes** completely protected from light at room temperature
- add **50µl of stop solution** in the same order as chromogen
- determine the optical extinction at **450nm** within 30 min after completing the assay.

Evaluation of the Test:

The extinctions determined for the calibration points are marked versus the concentrations of the standard solutions semilogarithmically (3.125 APL; 12.5 APL; 50 APL; 100 APL). A calibration curve is thus obtained by drawing a best-fit line (four parameters) through the points. Directly from the calibration curve, read the concentration of each analyte.

Cardiolipin-IgA values lower than 6 APL are judged as negative.

Cut-off = 6 APL

Interpretation of results depends on the specific clinical application of the test: each laboratory should establish its own clinically relevant ranges for the population taken into consideration. A specimen with equivocal antibody levels should be tested again; if it remains equivocal, the result should be reported as equivocal and/or an additional sample should be taken for testing according to physician judgment.

The test can be evaluated if the positive control serum is determined in the range given by the data sheet and if in the same time the negative control is found below the "Cut-Off" value.

Precautions:

For in-vitro diagnostics only! Standard and control sera are of human origin. The sera were tested and found to be negative for HBsAg, Hepatitis C, and HIV. Nevertheless all reagents should be regarded as potential infectious and thus handled with the care required. Rules for handling human sera have to be obeyed.

Warning: Some reagents contain Sodium Azide. Sodium Azide can form explosive Metal Azides with plumb and copper. Residuals of reagents should be removed with water carefully. Some reagents contain small amounts of Thimerosal (<0.1% w/v). Substrate contains 3,3', 5,5' Tetramethylbenzidine. The stop solution contains 2,6% Sulphuric Acid. Therefore handle all components as if potentially hazardous.

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- CRONIN ME, BISWAS RM, VAN DER STRAETON C, FLEISCHER TA et al.: IgG and IgM anticardiolipin antibodies in lupus anticardiolipin antibody associated clinical syndromes. *J. Rheumatol.*, 15: 795-798 (1988).
- COWCHOCK S, FORT J, MUNOZ S, NORBERG R et al. False positive ELISA test for anticardiolipin antibodies in sera from patients with repeated abortions, rheumatologic disorders, and primary biliary cirrhosis: Correlation with elevated polyclonal IgM and implications for patients with repeated abortions. *Clin. Exp. Immunol.*, 73: 289-294 (1988).
- HARRIS EN. Solid-phase anticardiolipin test revisited [editorial], *Am. J. Med.*, 85: 599-601 (1988).
- WILSON WA; GHARAVI AE, KOIKE T, LOCKSHIN MD et al. International consensus statement on preliminary classification criteria for definite antiphospholipid syndrome. *Arthritis & Rheumatism Vol.42 No.7:* 1309-1311 (1999)
- LOCKSHIN MD, SAMMARITANO LR, SCHWARTZMAN S. Validation of the Sapporo criteria for antiphospholipid syndrome. *Arthritis & Rheumatism Vol.43 No.2:* 440-443 (2000)
- MIYAKIS S, LOCKSHIN MD, ATSUMI T, BRANCH DW, BREY RL, CERVERA R et al. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J. Thromb and Haemost* 4: 295-306 (2006)

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- CRONIN ME, BISWAS RM, VAN DER STRAETON C, FLEISCHER TA et al.: IgG and IgM anticardiolipin antibodies in lupus anticardiolipin antibody associated clinical syndromes. *J. Rheumatol.*, 15: 795-798 (1988).
- COWCHOCK S, FORT J, MUNOZ S, NORBERG R et al. False positive ELISA test for anticardiolipin antibodies in sera from patients with repeated abortions, rheumatologic disorders, and primary biliary cirrhosis: Correlation with elevated polyclonal IgM and implications for patients with repeated abortions. *Clin. Exp. Immunol.*, 73: 289-294 (1988).
- HARRIS EN. Solid-phase anticardiolipin test revisited [editorial], *Am. J. Med.*, 85: 599-601 (1988).
- WILSON WA, GHARAVI AE, KOIKE T, LOCKSHIN MD et al. International consensus statement on preliminary classification criteria for definite antiphospholipid syndrome. *Arthritis & Rheumatism Vol.42 No.7*: 1309-1311 (1999)
- LOCKSHIN MD, SAMMARITANO LR, SCHWARTZMAN S. Validation of the Sapporo criteria for antiphospholipid syndrome. *Arthritis & Rheumatism Vol.43 No.2*: 440-443 (2000)
- MIYAKIS S, LOCKSHIN MD, ATSUMI T, BRANCH DW, BREY RL, CERVERA R et al. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J. Thromb and Haemost* 4: 295-306 (2006)

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CE Dr. Fenning BioMed GmbH