
EURO-DIAGNOSTICA

Instruction

Wieslab[®] Capture MPO-ANCA

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Performance characteristics

Table 1. Clinical sensitivity and specificity. A total of 348 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total number	Negative <5 IU/mL	Equivocal 5-7 IU/mL	Positive >7 IU/mL
Blood donors:	120	120	0	0
WG:	54	46	2	6
MP:	79	32	1	46
GBM:	55	42	4	9
SLE:	14	14	0	0
RA:	14	13	0	1
Sjögren's syndrome	12	12	0	0

WG = Wegener's granulomatosis, MP = microscopic polyangiitis
 SLE = systemic lupus erythematosus RA = rheumatoid arthritis
 GBM = glomerular basement membrane

Clinical sensitivity (equivocal samples excluded from calculations)

WG = $6/52 = 11.5\%$ 95% CI = 0.44 – 23.4%
 MP = $46/78 = 59.0\%$ 95% CI = 47.3 – 70.0%

Clinical specificity (equivocal samples excluded from calculations)

GBM = $42/51 = 82.4\%$ 95% CI = 69.1 – 91.6%
 SLE = $14/14 = 100\%$ 95% CI = 76.8 - 100%
 RA = $13/14 = 92.9\%$ 95% CI = 66.1 – 99.8 %
 Sjögren's syndrome = $12/12 = 100\%$ 95% CI = 73.5 - 100%
 Donors = $120/120 = 100\%$ 95% CI = 97.0 - 100%

The 95% confidence interval (CI) was calculated using the exact method.

Table 2. Agreement between the Capture MPO-ANCA kit and an alternative capture MPO-ELISA.

A total of 348 frozen retrospective sera were assayed. The following table summarises the results.

Wieslab® Capture-MPO-ANCA KIT					
Alternative ELISA		Positive	Equivocal	Negative	Total
	Positive	58	5	2	65
	Equivocal	0	0	1	1
	Negative	4	2	276	282
Total	62	7	279	348	

Sera falling in the equivocal range were excluded from the following calculations

Positive Percent Agreement: $58/60 = 96.7\%$ 95% CI = 88.5– 99.6%
 Negative Percent Agreement: $276/280 = 98.6\%$ 95% CI = 89.5 – 99.6%
 Overall Percent Agreement: $334/340 = 98.2\%$ 95% CI = 91.4 – 99.4%

The 95% confidence interval (CI) was calculated using the exact method.

Table 3. Batch to batch variation was determined by testing six different samples in eight replicates on three different batches.

	1	2	3	4	5	6
Mean value (IU/mL)	8	6	9	33	44	58
SD	0.85	0.74	0.62	1.69	2.38	4.35
% CV	11	13	7	5	5	8

Table 4. Inter-assay precision was determined by testing six different samples in eight replicates at three separate occasions.

	1	2	3	4	5	6
Mean value (IU/mL)	8	6	10	31	44	58
SD	0.57	0.74	1.61	5.20	6.31	6.86
% CV	7	11	17	17	14	12

Table 5. Intra-assay precision was determined by testing six different samples in eight replicates at one occasion.

	1	2	3	4	5	6
Mean value (IU/mL)	7	8	11	28	36	52
SD	0.49	0.93	0.69	1.46	3.07	2.21
% CV	7	11	6	5	9	4

Table 6. Dilution recovery was determined by testing five serial dilutions for three different samples.

Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution Corrected % Recovery
1	1/100	187	187	100
	1/200	101	94	108
	1/400	54	47	116
	1/800	28	23	120
	1/1600	14	12	120
Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution Corrected % Recovery
2	1/100	87	87	100
	1/200	45	44	103
	1/400	23	22	106
	1/800	12	11	110
	1/1600	6	5	110
Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution Corrected % Recovery
3	1/100	81	81	100
	1/200	43	41	106
	1/400	21	20	104
	1/800	10	10	99
	1/1600	5	5	99

Since patient sera contain heterogeneous antibody populations some samples may exhibit non-linearity, especially at very high sample dilutions.

Troubleshooting

Problem:	Possible causes:	Corrective action:
One or more calibrators out of range.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Improper pipetting. 	<ol style="list-style-type: none"> 1. Test invalid. Repeat test. 2. Test invalid. Repeat test.
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagent. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water solution used to prepare. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

References:

1. **Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): Manual of biological markers of disease. Dordrecht, Kluwer Academic Publishers; 1993:A9,1-14.
2. **Rasmussen N, Sjolín C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). J Immunol Methods 1990;127:139-145.
3. **Wieslander J.** How are anti neutrophil cytoplasmic autoantibodies detected? Am J Kidney Dis 1991;18:154-158.
4. **Jenne DE, Tschopp J, Ludemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. Nature 1990;346:520.
5. **Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
6. **Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
7. **Jennette C, Falk R.** Small vessels vasculitis. N Engl J Med 1997, 337, 1512-1525.
8. **Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. Kidney Int 1995, 48, 844-850.
9. **Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The T, van der Hem GK et al.** Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. Lancet 1990; 336:709-711.
10. **Savigne J, Davis D, Falk DJ, Jennette JC, Wiik A.** ANCA and associated diseases; A review of clinical and laboratory features. Kidney Int 2000, 57, 846-862.
11. **Savigne J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. Am J Clin. Pathol 1999, 111, 507-513.
12. **Csernok E, Ahlquist D, Ullrich S, Gross W.L.** A critical evaluation of commercial immunoassays for antineutrophil cytoplasmic antibodies directed against proteinase 3 and myeloperoxidase in Wegener's granulomatosis and microscopic polyangiitis.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklarung der Symbole. La spiegazione di simboli. Forklaring til symboler. Symboler som brukes p etiketter. Forklaringar till symboler.

	<p>Batch number. Numero de lot. El numero de la serie. Lot Nummer. Il numero di lotto. Lot nummer. Lot-nummer. Lot nummer.</p>
	<p>Product number. Reference. El numero del producto. Produktzahl. Il numero di prodotto. Bestillingsnummer. Katalognummer. Produktnummer.</p>
	<p>Expiration date. Date d'expiration. Fecha de caducidad. Verfalldatum. La data di scadenza. Udløbsdato. Brukes innen. Utgangsdatum.</p>
	<p>Store at. Conserver . Almacenar a. Lagerung bei. La conservatione. Opbevaringstemperatur. Temperaturbegrensning. Forvaringstemperatur.</p>
	<p>Biological material. Materiel biologique. Material biologico. Biologisches Material. Il materiale biologico. Biologisk materiale. Biologisk risiko. Biologiskt material.</p>
	<p>See instruction for use. Voir notice d'emploi. Siehe Anweisung zum Gebrauch. Vea la instruccion para el uso. Vedere l'istruzione per l'uso. Se brugsanvisning. Les bruksanvisningen. Se bruksanvisning.</p>
	<p>In vitro diagnostic medical device. Dispositif medical pour diagnostic in vitro. En vitro artefacto medico diagnostico. In Vitro diagnostische medizinische Vorrichtung. In vitro diagnostico medico congegno. In vitro diagnostisk brug. Til In vitro diagnostisk bruk. In vitro diagnostika.</p>
	<p>Manufacturer. Fabricant. El fabricante. Hersteller. Il fabbricante. Producent. Tilvirker. Tillverkare.</p>
	<p>Number of tests. Nombre de tests. Numero de pruebas. Anzahl der Tests. Numero di test. Antall tester. Antal tester. Siffror i symbolen anger antalet test.</p>
<p>Xi, Xn</p> 	<p>Irritant, Harmful. Irritant, Nocif. Irritante, Nocivo. Reizend, Gesundheitsgefahrdend. Irritante, Nocivo. Irriterende, Helsefarlig. Irriterande, Halsofarlig. Irriterande, halsofarlig.</p>

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Fortynning. Diluent. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Washsolution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentret 30 ganger. Tvättbuffert 30x konc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
SOLN STOP	Stop solution. Solution d'Arrêt. Solución de Parada. Stopplösung. Soluzione bloccante. Stopløsning. Stoppløsning. Stopplösning
CAL X	Calibrator. Etalon. Calibrador. Calibratore. Kalibrator. Kalibrator. Kalibrator. Kalibrator.
CONTROL X	Control. Contrôle. Controllo. Kontrolle. Controllo. Kontroll. Kontroll. Kontroll.