**Instruction**

**Zenit a-GBM semi quantitative kit  
(Goodpasture syndrome)**

**Enzyme immunoassay for detection of autoantibodies  
against glomerular basement membrane**

Microtitration strips (12x8) 96 wells  
Store the kit at +2-8°C  
For in vitro diagnostic use only

Document No. E-23-0165-02 M  
June, 2009

**Zenit a-GBM semi quantitative kit  
(Goodpasture syndrome)**

English:	page .....	2
Français:	page .....	9
Español:	página .....	17
Deutsch:	Seite .....	21
Italiano:	pagina .....	25
Português:	pagina .....	29
Ελληνικά:	σελ. ....	33

**REF** 40890

**IVD**

96

## **INTENDED USE**

The Zenit anti GBM test kit is an enzyme-linked immunosorbent assay (ELISA) for detection and semi-quantitation of IgG antibodies to glomerular basement membrane (GBM) in human sera. The assay is used to detect antibodies in a single serum specimen. The results of the assay are to be used as an aid to the diagnosis of Goodpasture syndrome. The analysis should be performed by trained laboratory professionals.

FOR IN VITRO DIAGNOSTIC USE.

## **Summary and explanation**

Goodpasture syndrome is characterised by lung haemorrhage, renal failure and the presence of anti-GBM antibodies. Diagnosis is based on clinical signs of lung hemorrhage and rapidly progressive glomerulonephritis combined with findings of anti-GBM antibodies.

As there are several other autoimmune diseases which may present with similar symptoms, this kit is a useful aid in differentiation between these diseases.

Less than one third of patients with reno-pulmonary syndromes have antibodies against the Goodpasture antigen, the majority having either Proteinase 3-ANCA or Myeloperoxidase-ANCA (1). Historically, indirect immunofluorescence was used to detect anti-GBM antibodies. When the first ELISA based on a collagenase digest was published in 1981 (2), assays using crude extracts were the only alternative. In 1984 the specific antigen was shown to derive from the C-terminal domain of type IV collagen (alpha 3 chain) (3, 4, 5, 6, 7, 11, 12, 13) and sensitive and specific assays were subsequently developed (8). The molecular nature of the Goodpasture antigen is reviewed in (9, 14). The antigen is characterised by a restricted tissue distribution occurring mainly in kidney and lungs 10).

The sensitivity and specificity of ELISA to detect anti-GBM antibodies is high (15). False positive reactions occur, mainly in SLE and other diseases with polyclonal activation, most of which yield low titers of 10-20 U/ml.

## **Principle of the Zenit anti-GBM kit**

The wells of the microtitre strips are coated with purified GBM antigen (alfa 3 chain). During the first incubation, specific antibodies in diluted serum, will bind to the antigen coating.

The wells are then washed to remove unbound antibodies and other components.

A conjugate of alkaline phosphatase-labelled antibodies to human IgG binds to the antibodies in the wells in this second incubation.

After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour

intensity and is measured in terms of absorbance (optical density (OD)). The absorbance is then calculated against a calibrator curve and the results are given in arbitrary U/ml.

## **Warnings and precautions**

- For in vitro diagnostic use.
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain proclin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing proclin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of anti-GBM in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity

### **Specimen collection**

The anti-GBM assay is for use with serum. Handle as if capable of transmitting infectious agents. Avoid using sera which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8° C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20° C or colder. Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

### **Kit components and storage of reagents**

- One frame with strips (12x8) coated with GBM antigen (alfa 3 chain) one lid sealed in a foil pack with a dry pack.
- 1,5 ml negative control (NC) containing human serum in diluent.
- 1,5 ml positive control (PC) containing human serum in diluent.
- 13 ml conjugate containing alkaline phosphatase-labelled antibodies to human IgG (blue colour).
- 32 ml Diluent (Dil) containing PBS (red colour).
- 13 ml Substrate pNPP.
- 30 ml wash solution 30x concentrated.
- Five calibrators containing human serum in diluent. 1.5 ml Cal 1 = 320 U/ml, 1.5 ml Cal 2 = 160 U/ml, 1.5 ml Cal 3 = 80 U/ml, 1.5 ml Cal 4 = 40 U/ml, 1.5 ml Cal 5 = 10 U/ml.

All reagents in the kit are ready for use except wash solution and should be stored at 2-8°C.

Remove only the number of strips needed for testing, resealing the aluminium package carefully.

### **Materials or equipment required but not provided**

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.

### **PROCEDURE**

All solutions should be used at room temperature. Incubate in all steps at room temperature (20-25°C) with lid. Incubate serum for 30 minutes, conjugate for 30 minutes and substrate for 60 minutes ( $\pm$  10 minutes).

### **Preparation of washing solution**

Dilute 10 ml of the 30x concentrated wash solution in 290 ml distilled water. When stored at 2-8°C, the diluted wash solution is stable until the date of expiration of the kit.

### **Dilution of serum and incubation**

Dilute the patient's serum 1/80 with diluent (395  $\mu$ l diluent +5  $\mu$ l serum).

Pipet 100 $\mu$ l/well in duplicate of diluent (as a blank), Calibrator 1, 2, 3, 4, 5, NC, PC and diluted patient's serum (P) according to the diagram below. Incubate for 30 minutes.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Dil	Cal 4	P1									
B	Dil	Cal 4	P1									
C	Cal 1	Cal 5	P2									
D	Cal 1	Cal 5	P2									
E	Cal 2	NC	etc									
F	Cal 2	NC										
G	Cal 3	PC										
H	Cal 3	PC										

### After serum incubation

Wash 3 times with 300 µl washing solution / well, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

### Adding conjugate

Add 100 µl conjugate to each well. Incubate for 30 minutes.

### After conjugate incubation

Wash as before.

### Adding substrate solution

Add 100 µl substrate pNPP to each well, incubate for 60 minutes ( $\pm$  10 minutes).  
Read the absorbance at 405 nm on a microplate reader.


### Calculations

Subtract the OD value for the blank from the other OD values.

Construct a calibrator curve by plotting the OD against the U/ml values of the 5 calibrators. The five calibrators provided have been assigned arbitrary values of 320 U/ml for calibrator 1, 160 U/ml for calibrator 2, 80 U/ml for calibrator 3, 40 U/ml for calibrator 4 and 10 U/ml for calibrator 5. Read the U/ml value of the patient from the constructed curve. Values greater than 320 should be reported as >320, or reassay them with a higher dilution.

Arbitrary U/ml have been adopted by Zenit, as no generally recognised international standard exists for expressing GBM titres.

<u>Example:</u>	<u>Calibrator</u>	<u>U/ml</u>	<u>Absorbance</u>
	1	320	2.3
	2	160	1.8
	3	80	1.2
	4	40	0.7
	5	10	0.18

A sample with an absorbance value of 0.7 will read on the X-axis as having 40 U/ml of anti-GBM.

**Important:** The curve is an example and should not be used for actual patient interpretation.

### Quality Control

The OD for the negative control should be less than that of calibrator 5.

The OD for calibrator 1 should be > 1,0

The value for the positive and negative controls see lot certificate.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off. It is recommended to assay an additional control at the assay cut-off.

If any of the values are not within their respective ranges, the test should be considered invalid and the test should be repeated. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to NCCLS C24-A for guidance on appropriate QC practices.

### Interpretation of results

< 10 U/ml = **Negative**

10-20 U/ml = **Equivocal**; Retest, if still equivocal retest by an alternative method or test a new sample

>20 U/ml = **Positive**

### Limitations

The individual patient's antibody titre can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardisation of results.

The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used as an adjunct to clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive for anti-GBM antibodies with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies.

Immunosuppressive therapy should not be started on basis of a positive anti-GBM result. Initiation or changes in treatment should not be based on changes in anti-GBM concentration alone, but rather on careful clinical observation.

### Expected results

Anti-GBM are rarely found in normal healthy individuals. The anti-GBM ELISA was tested with 120 normal sera. 120 were found to be negative. The annual incidence of Goodpasture syndrome is one per million inhabitants per year.

Anti-GBM antibodies are found in the sera of most Goodpasture syndrome patients. The anti-GBM ELISA was tested with 40 GP sera, 38 (95%) were found to be positive.

### Performance characteristics

**Table 1. Clinical sensitivity and specificity.** A total of 199 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total number	Negative <10 U/ml	Equivocal 10-20 U/ml	Positive >20 U/ml
<b>Blood donors:</b>	80	80	0	0
<b>WG / MP</b>	27	27	0	0
<b>SLE:</b>	19	19	0	0
<b>Others:</b>	33	33	0	0
<b>Anti-GBM:</b>	40	0	2	38

WG = Wegener's granulomatosis, MP = microscopic polyangiitis

SLE = systemic lupus erythematosus Others = RA, UC, etc.

GBM = glomerular basement membrane

### Clinical sensitivity (Equivocal samples excluded from calculations)

GBM = 38/38 = 100 %

95% CI = 90.7-100 %

**Clinical specificity (Equivocal samples excluded from calculations)**

WG/MP	= 27/27 = 100 %	95% CI = 87.2-100 %
SLE	= 19/19 = 100 %	95% CI = 82.4-100 %
Others	= 33/33 = 100 %	95% CI = 89.4-100 %
Donors	= 80/80 = 100 %	95% CI = 95.5-100 %

The 95% confidence interval (CI) was calculated using the exact method.

**Table 2. Relative sensitivity and specificity of the Zenit anti-GBM kit compared to GBM IFA.**

A total of 68 frozen retrospective sera were assayed. The following table summarises the results.

<b>Zenit anti-GBM</b>				
<b>GBM IFA</b>	<b>Positive</b>	<b>Equivocal</b>	<b>Negative</b>	<b>Total</b>
Positive	55	3	0	58
Negative	1	1	8	10
<b>Total</b>	<b>56</b>	<b>4</b>	<b>8</b>	<b>68</b>

Sera falling in the equivocal range were excluded from the following calculations.

Relative sensitivity	= 55/55 = 100 %	95% CI = 93.5-100 %
Relative specificity	= 8/9 = 88.9 %	95% CI = 51.8-99.7 %
Relative accuracy	= 63/64 = 98.4%	95% CI = 91.6-100 %

**Table 3. Relative sensitivity and specificity of the Zenit anti-GBM kit compared to an alternative ELISA.**

(Zenit 60x60x60) A total of 159 frozen retrospective sera were assayed. The following table summarises the results.

<b>Zenit anti-GBM</b>					
		<b>Positive</b>	<b>Equivocal</b>	<b>Negative</b>	<b>Total</b>
<b>Alternative ELISA</b>	<b>Positive</b>	35	0	0	35
	<b>Equivocal</b>	3	1	0	4
	<b>Negative</b>	0	1	119	120
	<b>Total</b>	<b>38</b>	<b>2</b>	<b>119</b>	<b>159</b>

Sera falling in the equivocal range were excluded from the following calculations.

Relative sensitivity	= 35/35 = 100 %	95% CI = 90.0-100 %
Relative specificity	= 119/120 = 99.2 %	95% CI = 95.4-100 %
Relative accuracy	= 154/155 = 99.4%	95% CI = 96.5-100 %

The 95% confidence interval (CI) was calculated using the exact method.

**Table 4. Batch to batch variation** was determined by testing four different samples in duplicate. Results were obtained for four different batches.

<b>Sample</b>	<b>Mean value</b>	<b>SD</b>	<b>CV %</b>
1	58 U/ml	5	9
2	45 U/ml	7	15
3	24 U/ml	2	9
PK	40 U/ml	5	11

**Table 5. Inter-assay precision** was determined by testing two different samples in duplicate. Results were obtained for seven different runs.

Sample	Mean value	SD	CV %
1	40 U/ml	2.6	7
2	18 U/ml	1.9	11

**Table 6. Intra-assay precision** was determined by testing one samples in 70 wells.

Sample	Mean value	SD	CV %
1	67 U/ml	5.4	8

**Table 7. Linearity.** The values were determined for serial two-fold dilutions of four positive sera. The values were compared to log 2 of dilution by standard linear regression. The data indicates that the assay has a linear relationship with serum dilution.

Serum	neat	1:2	1:4	1:8	1:16	r
1	60	36	18	9		0,977
2	217	126	82	71	26	0,988
3	98	60	34	23	6	0,998
4	69	42	19	12		0.999

## Troubleshooting

<b>Problem:</b>	<b>Possible causes:</b>	<b>Solution:</b>
Control values out of range.	<ol style="list-style-type: none"> <li>1. Incorrect temperature, timing or pipetting; reagents not mixed.</li> <li>2. Cross contamination of controls.</li> <li>3. Improper dilution.</li> <li>4. Optical pathway not clean.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test.</li> <li>2. Pipette carefully.</li> <li>3. Repeat test.</li> <li>4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.</li> </ol>
All test results negative.	<ol style="list-style-type: none"> <li>1. One or more reagents not added, or added in wrong sequence.</li> <li>2. Antigen coated plate inactive.</li> </ol>	<ol style="list-style-type: none"> <li>1. Recheck procedure. Check for unused reagents. Repeat test.</li> <li>2. Check for obvious moisture in unused wells. Wipe bottom and reread.</li> </ol>
All test results yellow.	<ol style="list-style-type: none"> <li>1. Contaminated buffers or reagents.</li> <li>2. Washing solution contaminated</li> <li>3. Improper dilution of serum.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check all solutions for turbidity.</li> <li>2. Use clean container. Check quality of water used to prepare solution.</li> <li>3. Repeat test.</li> </ol>
Poor precision.	<ol style="list-style-type: none"> <li>1. Pipette delivery CV greater than 5%</li> <li>2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature.</li> <li>3. Reagent addition taking too long; inconsistency in timing intervals.</li> <li>4. Optical pathway not clean.</li> <li>5. Washing not consistent; trapped bubbles; washing solution left in the wells</li> <li>6. Improper pipetting.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check calibration of pipette. Use reproducible technique.</li> <li>2. Mix all reagents gently but thoroughly and equilibrate to room temperature.</li> <li>3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time.</li> <li>4. Check for air bubbles in the wells. Wipe bottom and reread.</li> <li>5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue.</li> <li>6. Avoid air bubbles in pipette tips.</li> </ol>

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**REFERENCES. RÉFÉRENCES. REFERENCIAS. LITERATUR. BIBLIOGRAFIA. ΑΝΑΦΟΡΕΣ:**

1. **Saxena R, Bygren P, Arvastson B, Wieslander J.** Circulating autoantibodies as serological markers in the differential diagnosis of pulmonary renal syndrome. *J Intern Med.* 1995. 238: 143-152.
2. **Wieslander J, Bygren P, Heinegård D.** Anti-basement membrane antibody: Immunoenzymatic assay and specificity of antibodies. *Scand. J. Clin. Lab. Invest.* 1981. 41:763-772.
3. **Wieslander J, Bygren P, Heinegård D.** Isolation of the specific glomerular basement membrane antigen involved in Goodpasture syndrome. *Proc Natl Acad Sci U. S. A.* 1984a. 81:1544-1548.
4. **Wieslander J, et al** Goodpasture antigen of the glomerular basement membrane: localization to noncollagenous regions of type IV collagen. *Proc Natl Acad Sci U. S. A.* 1984b. 81:3838-3842.
5. **Segelmark M, Butkowski R, Wieslander J.** Antigen restriction and IgG subclasses among anti-GBM autoantibodies. *Nephrol Dial Transplant* 1990.5: 991-996.
6. **Hellmark T, Johansson C, Wieslander J.** Characterization of anti-GBM antibodies involved in Goodpastures syndrome. *Kidney Int.* 1994. 46: 823-829.
7. **Hellmark T, Brunmark C, Trojner J, Wieslander J.** Epitope mapping of anti-GBM antibodies with synthetic peptides. *Clin. Exp. Immunol.* 1996, 105 504-510.
8. **Saxena R, Isaksson B, Bygren P, Wieslander J.** A rapid assay for circulating anti-glomerular basement membrane antibodies in Goodpasture syndrome. *J. Immunol. Methods* 1989 118:73-78.
9. **Hellmark T, Segelmark M, Bygren P, Wieslander J.** Glomerular basement membrane antibodies. In "Autoantibodies" Eds Peter JB, Shoenfeld Y, Elsevier Science 1996, 291-298.
10. **Kleppel MM, et al.** Human tissue distribution of novel basement membrane collagen. *Am J Path.* 1989. 134 4: 813-825.
11. **Hellmark T, et al.** Identification of a clinically relevant immunodominant region in Goodpasture disease. *Kidney Int* 1999, 55, 936-944.
12. **Segelmark M, Burkhardt H, Wieslander J.** Goodpasture disease: Characterization of a single conformational epitope as the target of pathogenic autoantibodies. *J. Biol. Chem.* 1999, 274, 25862-25868.
13. **Gunnarsson A, et al.** Molecular properties of the Goodpasture epitope. *J Biol Chem* 2000, 275, 30844-30848.
14. **Hudson B, et al..** Alport's Syndrome Goodpasture's Syndrome and type IV collagen *N Engl J Med* 2003, 348, 2543-2556.
15. **Segelmark M, Hellmark T, Wieslander J..** The prognostic significance in Goodpasture's disease of specificity, titre and affinity of anti-glomerular basement membrane antibodies. *Nephron Clin Pract* 2003, 94, 59-68.

**EXPLANATION OF SYMBOLS. EXPLICATION DE SYMBOLES. EXPLICACIÓN DE SIMBOLOS. ERKLÄRUNG DER SYMBOLE. SPIEGAZIONE DEI SIMBOLI. EXPLICAÇÃO DOS SÍMBOLOS. ΕΠΕΞΗΓΗΣΗ ΣΥΜΒΟΛΩΝ.**

	Batch number. Numéro de lot. Número de lote. Chargennummer. Numero di lotto. Código do lote. Αριθμός παρτίδας.
	Product number. Référence. Número del producto. Produktnummer. Numero di prodotto. Referência de catálogo. Αριθμός προϊόντος
	Expiration date. Date d'expiration. Fecha de caducidad. Verfallsdatum. Data di scadenza. Prazo de validade. Ημερομηνία λήξης.
	Store at. Conserver à. Conservar a. Lagerung bei. Conservare a. Conservar a. Φυλάσσετε σε:
	Biological material. Matériel biologique. Material biológico. Biologisches Material. Materiale biologico. Material biológico. Βιολογικό υλικό.
	See instruction for use. Voir notice d'emploi. Consultar las instrucciones de uso. Gebrauchsanweisung beachten. Consultare le istruzioni per l'uso. Consulte as instruções de utilização. Βλ. Οδηγίες χρήσης.
	In vitro diagnostic medical device. Dispositif médical pour diagnostic in vitro. Dispositivo médico para diagnóstico <i>in vitro</i> . In-vitro-Diagnostikum. Dispositivo medico diagnostico in vitro. Dispositivo médico para diagnóstico <i>in vitro</i> . Ιατροτεχνολογικό βοήθημα που χρησιμοποιείται στη διάγνωση in vitro.
	Manufacturer. Fabricant. Fabricante. Hersteller. Fabbricante. Fabricante. Κατασκευαστής
	Number of tests. Nombre de tests. Número de pruebas. Anzahl der Tests. Numero di test. Número de ensaios. Αριθμός εξετάσεων
	Irritant. Irritant. Irritante. Reizend. Irritante. Irritante. Ερεθιστικό.

<b>Ag</b>	Antigen. Antigène. Antígeno. Antigen. Antigene. Antigénio. Αντιγόνο.
<b>DIL</b>	Diluent. Diluant. Diluyente. Verdünnungsmittel. Diluente. Diluente. Διαλύτης
<b>CONJ</b>	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Συζυγές.
<b>BUF</b> <b>WASH</b> <b>30X</b>	Washsolution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konz. Soluzione di lavaggio conc. 30x. Solução de lavagem concentrada 30x. Διάλυμα πλύσης με συμπ. 30x.
<b>SUBS</b> <b>pNPP</b>	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Υπόστρωμα pNPP.
<b>SOLN</b> <b>STOP</b>	Stop solution. Solution d'Arrêt. Solución de interrupción. Stopplösung. Soluzione bloccante. Solução de paragem. Ανασχετικό διάλυμα.
<b>CAL</b> <b>X</b>	Calibrator. Etalon. Calibrador. Kalibrator. Calibratore. Calibrador. Βαθμονομητής.
<b>CONTROL</b> <b>X</b>	Control. Contrôle. Control. Kontrolle. Controllo. Controllo. Μάρτυρας

**UK**  
**UNITED KINGDOM**  
*Distributed by*  
A. Menarini Diagnostics Ltd  
405 Wharfedale Road  
Winnersh - Wokingham  
Berkshire RG41 5RA

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A. Menarini Diagnostics S.A.
Benelux S.A./N.V.
Chemie AG
Euro-Diagnostica AB, Malmo, Sweden