




unyvero

Unyvero's sample-to-answer platform provides rapid results for severe infectious diseases in hospitalized patients

Powerful multiplex PCR technology combined with the broadest range of microorganism and resistance targets sets the Unyvero System apart.


The Unyvero System consists of:

- Lysator to lyse and process a variety of native samples
- Cockpit to manage testing process, display, store, and transmit results
- Analyzer to perform DNA testing with random-access, multiplex PCR



A single test handles one patient sample, analyzes over 100 DNA analytes and delivers reliable results within just 4-5 hours

unyvero
a product by curetis




Unyvero is designed to expand with your growing needs

Applications for severe infections:

- Blood Culture – BCU
- Hospitalized Pneumonia – HPN
- Intra-Abdominal Infection – IAI
- Implant & Tissue Infection – ITI
- Urinary Tract Infection – UTI



unyvero
a product by curetis



Implant & Tissue Infection

Fast & Simple Syndromic Testing for Severe Infections - Improving Patient Outcomes



Unyvero L4 Lysator



Unyvero C8 Cockpit



Unyvero A50 Analyzer

The Unyvero System is distributed on an exclusive basis by A.Menarini Diagnostics in the following countries: Benelux, France, Germany, Greece, Italy, Portugal, Spain, United Kingdom.

Distributed by: A.Menarini Diagnostics S.r.l.

Legal Site: Via Sette Santi, 3 - 50131 Florence, Italy - Operative Site: Via Lungo L'Emma, 7 - 50012 Bagno a Ripoli (FI), Italy - Website: www.menarini.com

- BENELUX - MENARINI BENELUX N.V./S.A: Division Diagnostics De Kleetlaan 3/1831 Machelen/Belgium
Tel. 0032-2-721 49 30 - Fax 0032-2-721 50 49
- FRANCE - A. MENARINI DIAGNOSTICS FRANCE S.A.S.U.: 3 - 5, rue du Jura BP 70511 F-94633 Rungis Cédex
Tel. 0033-1-56 34 69 10 - Fax 0033-1-56 34 69 11
- GERMANY - A. MENARINI DIAGNOSTICS DEUTSCHLAND: Division der BERLIN-CHEMIE AG Glienicke Weg 125 D-12489 Berlin
Tel. 0049-30-6707 3000 Fax 0049-30-6707 3020
- GREECE - MENARINI DIAGNOSTICS S.A.: 575, Vouliagmenis Ave. GR-164 51 Argypolis - Athens
Tel. 0030-210-99 44 952 - Fax 0030-210-99 45 029
- ITALY - A. MENARINI DIAGNOSTICS S.r.l.: Via Lungo l'Emma, 7 50012 Bagno a Ripoli (Firenze)
Tel. 0039-055-5680-422/394 - Fax 0039-055-5680-905
- PORTUGAL - MENARINI DIAGNOSTICOS LDA: Quinta da Fonte Edifício D, Manuel I, 2º B 2770-203 Paço de Arcos - Lisbon
Tel. 00351-210-930 000 - Fax 00351-210-930 001
- SPAIN - MENARINI DIAGNOSTICOS S.A.: Avenida del Maresme, 120 08918 Badalona (Barcelona)
Tel. 0034-93-50 71 000 - Fax 0034-93-27 80 215
- UK - A. MENARINI DIAGNOSTICS LTD: 405 Wharfedale Road Wetherby, Wokingham Berkshire RG41 5RA
Tel. 0044-118 944 4100 - Fax 0044-118 944 4111

Caution - Investigational Device. Limited by Federal US Law to investigational Use. Not available for Sale in the United States.

FOR PROFESSIONAL USE ONLY Cod.51259



Implant and tissue infections can be complex and time-consuming to diagnose

- Treatment of patients with implant and tissue infection can only be optimized after the causative microorganism and its associated resistance are known.
- Each sample collected is analyzed through conventional microbiology. Some of them are cultured for more than 14 days.¹

In the UK, average cost of a knee revision for infection is estimated to be around £30,000.²

- Empiric broad spectrum antibiotics may not provide optimal coverage and can exacerbate resistance.³
- Biofilm formation can often develop on orthopedic implants and is difficult to diagnose with culture methods.⁴

¹ Bossard DA et al., Optimal Length of Cultivation Time for Isolation of Propionibacterium acnes in Suspected Bone and Joint Infections Is More than 7 Days. J Clin Microbiol. 2016;54(12):3043-3049.
² Kallala RF et al., Financial analysis of revision knee surgery based on NHS tariffs and hospital costs: does it pay to provide a revision service? Bone Joint J. 2015;97-B(2):197-201.
³ WHO Antimicrobial Resistance. Global Report on Surveillance. 2014.
⁴ Torres C et al., ECCMID 2017.

Faster detection enables earlier optimization of therapy

The Unyvero ITI Application includes a broad spectrum of microorganisms and resistance genes.

Unyvero ITI can be used for the diagnosis of:

- Burn wound infections
- Cardiology-associated infections
- Catheter-associated infections
- Deep skin and tissue infections
- Diabetic foot infections
- Orthopedic implant infections
- Surgical site infections

Clinical evidence demonstrates the benefits provided by the Unyvero solution

Study 1

Comparison between Unyvero ITI and conventional methods focusing on cost-analysis. Hospital Clínic de Barcelona.


Eligibility

Patients who underwent implant removal due to suspicion of infection.

Study design:

- Sonication of removed implants (hip 46%, knee 42%, shoulder 12%).
- Comparison of patient management cost between traditional microbiology (n=10) vs Unyvero (n=14).
- Cost analysis model based on antibiotic treatment (empiric and specific), hospital stay and Unyvero costs.

Torres et al., poster presented at European Congress of Clinical Microbiology (ECCMID) 2017.



Conclusion

Unyvero ITI supports a rapid diagnosis of PJI when an infection is suspected. Its use is associated with a shorter hospital length of stay compared to standard culture methods allowing cost savings at hospital level.

Study 2

Clinical evaluation of the ITI Application in patients suspected of prosthetic joint infection. Helios Endo Klinik Hamburg.

Number of samples

60 intra-operative joint aspirates.

Eligibility

Patients >18 years old undergoing hip or knee revision arthroplasty.

88.3% Accuracy

100% Specificity

Clinical performance

- Aseptic loosening (n=26):
Concordance = 100%.
- Chronic PJI (n=26):
Sensitivity = 76.9%,
Specificity = 100%.
- Acute PJI (n=8):
Sensitivity = 85.7%,
Specificity = 100%.

Time to results


- Mean time for conventional culture results was 6.4 days (range 48-552 hours) whereas Unyvero results were available in 5 hours.

148 hours saved

Systemic Inflammatory Response Syndrome

- 2 patients with SIRS.
- Joint aspirate culture negative vs positive with Unyvero ITI (1 CoNS, 1 E. cloacae).
- ITI results later confirmed with tissue culture.

Lausmann et al., 2017 J Bone Jt Infect 2(4): 175-183.



Unyvero Implant & Tissue Infection (ITI) Cartridge

Gram-positive bacteria	Enterobacteriaceae	Non-fermenting bacteria	Corynebacteriaceae	Resistance	Gene
<i>Staphylococcus aureus</i> Coagulase negative staphylococci <i>Streptococcus</i> spp. <i>Streptococcus agalactiae</i> <i>Streptococcus pneumoniae</i> <i>Streptococcus pyogenes/dysgalactiae</i> <i>Granulicatella adiacens</i> <i>Abiotrophia defectiva</i> <i>Enterococcus</i> spp. <i>Enterococcus faecalis</i>	<i>Citrobacter freundii/koseri</i> <i>Escherichia coli</i> <i>Enterobacter cloacae</i> complex <i>Klebsiella aerogenes</i> (E. aerogenes) <i>Klebsiella pneumoniae</i> <i>Klebsiella oxytoca</i> <i>Klebsiella varicola</i> <i>Proteus</i> spp. Universal bacteria Detection of prokaryotic genetic sequence	<i>Acinetobacter baumannii</i> complex <i>Pseudomonas aeruginosa</i> Anaerobic bacteria <i>Cutibacterium acnes</i> (P. acnes) <i>Fingoldia magna</i> <i>Bacteroides fragilis</i> group	<i>Corynebacterium</i> spp. Fungi <i>Candida</i> spp. <i>Candida albicans</i> <i>Candida glabrata</i> <i>I. orientalis</i> (C. krusei) <i>Candida tropicalis</i>	Macrolide/ Lincosamide Aminoglycoside Oxacillin Vancomycin 3rd generation Cephalosporins Carbapenem	<i>ermA</i> <i>ermC</i> <i>aac(6)/aph(2')</i> <i>aacA4</i> <i>mecA</i> <i>mecC</i> <i>vanA</i> <i>vanB</i> <i>ctx-M</i> <i>kpc</i> <i>imp</i> <i>ndm</i> <i>oxa-23</i> <i>oxa-24/40</i> <i>oxa-48</i> <i>oxa-58</i> <i>vim</i>

Sample Types

Sonication fluids, swabs, tissue, pus, aspirate/exudate, bone fragments, etc.



Easy Workflow



Multiple Sample Types



24/7 Results