




unyvero

Unyvero's sample-to-answer platform provides rapid results for severe infectious diseases in hospitalized patients

Powerful multiplex PCR technology combined with the broadest range of microorganism and resistance targets sets the Unyvero System apart.

The Unyvero System consists of:

- Lysator to lyse and process a variety of native samples
- Cockpit to manage testing process, display, store, and transmit results
- Analyzer to perform DNA testing with random-access, multiplex PCR



A single test handles one patient sample, analyzes over 100 DNA analytes and delivers reliable results within just 4-5 hours



Unyvero L4 Lysator



Unyvero C8 Cockpit



Unyvero A50 Analyzer

unyvero

a product by curetis



Unyvero is designed to expand with your growing needs

Applications for severe infections:

- Blood Culture – BCU
- Hospitalized Pneumonia – HPN
- Intra-Abdominal Infection – IAI
- Implant & Tissue Infection – ITI
- Urinary Tract Infection – UTI



The Unyvero System is distributed on an exclusive basis by A.Menarini Diagnostics in the following countries: Benelux, France, Germany, Greece, Italy, Portugal, Spain, United Kingdom.

Distributed by: A.Menarini Diagnostics S.r.l.

Legal Site: Via Sette Santi, 3 - 50131 Florence, Italy - Operative Site: Via Lungo L'Emma, 7 - 50012 Bagno a Ripoli (FI), Italy - Website: www.menarini.com


- BENELUX - MENARINI BENELUX N.V./S.A: Division Diagnostics De Kleetlaan 3/1831 Machelen/Belgium Tel. 0032-2-721 49 30 - Fax 0032-2-721 50 49
- FRANCE - A. MENARINI DIAGNOSTICS FRANCE S.A.S.U.: 3 - 5, rue du Jura BP 70511 F-94633 Rungis Cédex Tel. 0033-1-56 34 69 10 - Fax 0033-1-56 34 69 11
- GERMANY - A. MENARINI DIAGNOSTICS DEUTSCHLAND: Division der BERLIN-CHEMIE AG Glienicke Weg 125 D-12489 Berlin Tel. 0049-30-6707 3000 Fax 0049-30-6707 3020
- GREECE - MENARINI DIAGNOSTICS S.A.: 575, Vouliagmenis Ave. GR-164 51 Argypoulis - Athens Tel. 0030-210-99 44 952 - Fax 0030-210-99 45 029
- ITALY - A. MENARINI DIAGNOSTICS S.r.l.: Via Lungo l'Emma, 7 50012 Bagno a Ripoli (Firenze) Tel. 0039-055-5680-422/394 - Fax 0039-055-5680-905
- PORTUGAL - MENARINI DIAGNOSTICOS LDA: Quinta da Fonte Edifício D, Manuel I, 2º B 2770-203 Paço de Arcos - Lisbon Tel. 00351-210-930 000 - Fax 00351-210-930 001
- SPAIN - MENARINI DIAGNOSTICOS S.A.: Avenida del Maresme, 120 08918 Badalona (Barcelona) Tel. 0034-93-50 71 000 - Fax 0034-93-27 80 215
- UK - A. MENARINI DIAGNOSTICS LTD: 405 Wharfedale Road Wetherby, Wokingham Berkshire RG41 5RA Tel. 0044-118 944 4100 - Fax 0044-118 944 4111

Caution – Investigational Device. Limited by Federal US Law to investigational Use. Not available for Sale in the United States.

FOR PROFESSIONAL USE ONLY Cod. 51258

unyvero

a product by curetis



# Intra-Abdominal Infection

Fast & Simple Syndromic Testing for Severe Infections - Improving Patient Outcomes



## Intra-abdominal infections (IAI) are often associated with poor prognosis, particularly in high-risk patients

- Early clinical diagnosis and appropriate antimicrobial therapy are essential in the management of intra-abdominal infections.<sup>1</sup>
- Conventional microbiology methods can take 1-4 days and diagnosis of anaerobic bacteria up to 14 days.
- Delayed or inadequate antimicrobial therapy is associated with poorer outcomes and increased death.<sup>2</sup>

IAI is the second most commonly identified cause of sepsis in the intensive care unit.<sup>3,4</sup>

- Global antibiotic resistance continues to rise steadily and can make therapy selection difficult.<sup>5</sup>
- Cost of antifungal treatment is very high.<sup>6</sup>



## Faster detection enables earlier optimization of therapy

The Unyvero IAI Application simultaneously detects a large panel of bacteria, fungi, toxins and antibiotic resistance markers directly from IAI samples.

The Unyvero IAI Cartridge can be used for the diagnosis of:

- Acute abdomen
- Ascites
- Cholecystitis
- Diverticulitis
- Peritonitis
- Surgical site infections



## Clinical evidence demonstrates the benefits provided by the Unyvero solution

### Study 1

CE Performance evaluation study.

Number of samples  
332 culture positive samples.

Sample types  
Ascites, bile, drainage fluid, pancreatic juice, peritoneal fluid, pus, swabs, tissue; positive blood cultures inoculated with ascites.

Additional identification

- Additional potential pathogens were detected in 308 cases.
- >292 were confirmed as true positives (PCR and sequencing).

Unyvero turnaround time 4-5 hours

Analytical performance

- All analytes included in the panel achieved a sensitivity of at least 75% while 13/27 targets achieved a sensitivity of 100%.
- Specificity ranged from 98.5% to 100% over all panel targets.

89.2%  
Sensitivity

99.5%  
Specificity

### Study 2

Multicentre Clinical Evaluation. Barts Health NHS Trust, The Great Romagna Hub, CHU Amiens and CHU Toulouse.

Number of samples  
300 samples (mainly peritoneal fluids, ascites, pus and bile).

Eligibility  
Samples from patients with suspicion of intra-abdominal infection.

Improved microorganism detection

- Additional microorganisms identified, in particular anaerobes, with most detections (91.4%) confirmed by sequencing.
- Microbiology reported results negative for 86 samples:

Of these, 62 samples (72%) also negative with Unyvero IAI while pathogens were detected in the remaining 24.

In 16/24 samples, pathogens detected by Unyvero IAI were confirmed by sequencing.

Time saving

- Using Unyvero, the average time to results was reduced by 17h compared to identification results.

17 hours saved



(39:06 ± 16:09h for microbiology vs 22:02 ± 4:12h Unyvero IAI).

- Using Unyvero, the average time to results was reduced by 41h compared to full AST results.

41 hours saved



(64:19 ± 12:10h for microbiology vs. 23:44 ± 3:58h Unyvero IAI).

## Unyvero Intra-Abdominal Infection (IAI) Cartridge

Gram-positive bacteria	Enterobacteriaceae	Anaerobic/facultative anaerobic bacteria	Non-fermenting bacteria	Resistance	Gene
<i>Coagulase negative staphylococci</i> <i>Enterococcus faecalis</i> <i>Enterococcus</i> spp. <i>Streptococcus</i> spp. <i>Staphylococcus aureus</i>	<i>Escherichia coli</i> <i>Klebsiella aerogenes</i> ( <i>E. aerogenes</i> ) <i>Enterobacter cloacae</i> complex <i>Klebsiella pneumoniae</i> <i>Klebsiella oxytoca</i> <i>Klebsiella variicola</i> <i>Proteus</i> spp.	<i>Aeromonas</i> spp. <i>Bacteroides fragilis</i> group <i>Bacteroides</i> spp. / <i>Prevotella</i> spp. <i>Clostridium difficile</i> <i>Clostridium perfringens</i> <i>Fingoldia magna</i> <i>Cutibacterium acnes</i> ( <i>P. acnes</i> )	<i>Acinetobacter baumannii</i> complex <i>Pseudomonas aeruginosa</i>	Oxacillin Vancomycin Aminoglycoside 3rd generation Cephalosporins Fosfomycin Polypeptides/polymyxins Nitroimidazole Fluoroquinolone Tetracycline	<i>mecA</i> <i>mecC</i> <i>vanA</i> <i>vanB</i> <i>aac A4</i> <i>ctx-M</i> <i>fosA3</i> <i>mcr-1</i> <i>nimA</i> <i>nimB</i> <i>qnrA</i> <i>qnrB</i> <i>qnrS</i> <i>tetA</i>
	Universal bacteria	Toxin	Fungi		<i>kpc</i> <i>imp</i> <i>ndm</i> <i>oxa-23</i> <i>oxa-24/40</i> <i>oxa-48</i> <i>oxa-58</i> <i>vim</i>
	Detection of prokaryotic genetic sequence	Toxin B Shiga Toxin	<i>Candida</i> spp. <i>Candida albicans</i> <i>Candida glabrata</i> <i>Candida tropicalis</i> <i>Issatchenkia orientalis</i> ( <i>Candida krusei</i> )		
		Marker <i>tcdB</i> <i>stx1/2</i>			

### Sample Types

Ascites and peritoneal fluid, pancreatic juice, bile, tissue, puncture fluid, swabs, catheter/drainage tips, and samples from positive blood culture bottles that have been inoculated with ascites/puncture fluid.



Easy  
Workflow



Multiple  
Sample Types



24/7  
Results

<sup>1</sup> Sartelli M et al., Management of intra-abdominal infections: recommendations by the WSES 2016 consensus conference. World J Emerg Surg. 2017; 12:22.  
<sup>2</sup> Solomkin JS et al., Diagnosis and management of complicated intra-abdominal infection in adults and children: guidelines by the Surgical Infection Society and the Infectious Diseases Society of America. Clin Infect Dis. 2010;50(2):133-64.  
<sup>3</sup> Lopez N et al., A Comprehensive review of abdominal infections. World J Emerg Surg. 2011; 6:7.  
<sup>4</sup> Sartelli Met al., Current concept of abdominal sepsis: WSES position paper. World J Emerg Surg. 2014; 9:22.  
<sup>5</sup> WHO Antibiotic Resistance Fact Sheet. November 2017.  
<sup>6</sup> Pagès A et al., Cost Effectiveness of Candida Polymerase Chain Reaction Detection and Empirical Antifungal Treatment among Patients with Suspected Fungal Peritonitis in the Intensive Care Unit. Value Health. 2017;20(10):1319-1328.

Ciesielczuk et al., Multicenter performance evaluation of the Unyvero IAI cartridge for detection of intra-abdominal infections. Eur J Clin Microbiol Infect Dis. 2018; 37(11):2107-2115